
anti- dog IgG antibody

Product Information

Catalog No.:	FNab10899
Size:	1mg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

This rabbit polyclonal antibody recognizes the dog IgG.

Immunogen information

Immunogen:	dog IgG
Synonyms:	dog IgG,canine IgG

Application

Reactivity:	dog
Tested Application:	ELISA, WB

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com