

anti- ASCT2 antibody

Product Information

Catalog No.:	FNab10940
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Solute carrier family 1 member 5
Synonyms:	AAAT, ASCT2, ATB(0), ATBO, Baboon M7 virus receptor, M7V1, M7VS1, R16, RDR, RDRC, SLC1A5
Observed MW:	70 kDa
Uniprot ID :	Q15758

Application

Reactivity:	Human, Mouse
Tested Application:	WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50 - 1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human colon cancer using FNab10940(ASCT2 Antibody) at dilution of 1:100



Various lysates were subjected to SDS PAGE followed by western blot with FNab10940(ASCT2 Antibody) at dilution of 1:1000.