



## Wuhan Fine Biotech Co., Ltd




+86-027-87384275 +86-027-87800889

fine@fn-test.com www.fn-test.com

B9 Bld, High-tech Medical Devices Park, No. 818 Gaoxin Ave.  
East Lake High-tech Development Zone  
Wuhan, Hubei, China(430206)

# YOUR GOOD RESEARCH PARTNER

Wuhan Fine Biotech Co., Ltd




## About Us



Wuhan Fine Biotech Co., Ltd.(FineTest®) is an ISO 9001 certificated manufacturer which focus on supplying bio-reagents for research use. We are dedicated to providing 7000+ types ELISA Kits, 10000+ primary & secondary antibodies, 4500+ recombinant proteins and other related reagents like assay kits to global customers in the research fields of immunology, cell biology, tumor, neurobiology etc.

With stable product quality, professional tech support and good service, our products are trusted by our faithful partners and customers around the world. Thousands of FineTest research products are cited by articles published in famous journals, such as Nature, Cell, Hepatology, Brain...etc.

## Product Research Areas



## ELISA Kits

The ready-to-use ELISA kits are our most trusted product line. There are 7000+ types involving in about 16 species not only regular species like Human, Rat, Mouse, but also rare species like Bovine, Chicken, Canine...etc., which can meet various kinds of research subject requirement.

**Ready-to-Use:** Include all reagents and components required by the ELISA test

- Pre-coated ELISA Microplate
- Biotin-Labeled Antibody
- HRP-Streptavidin Conjugate
- Wash Buffer
- Standard
- Optimized Sample Dilution Buffer
- TMB Substrate
- Stop Solution
- Plate Sealer
- Product Description

**Precision:** Ensured by validated recovery and linearity.

Linearity data of EH0201 Human IL-6 ELISA Kit.

Serum samples	Dilution	Average of Expected	Range
N=10	1:2	95-110%	99%
	1:4	91-107%	95%
	1:8	88-115%	92%
Cell culture samples	Dilution	Average of Expected	Range
N=10	1:2	97-120%	100%
	1:4	92-117%	97%
	1:8	90-119%	93%

**High Specificity:** Select high-specific antibody pairs and testify cross reaction between target and similar protein.

**Reproducibility:** Strict QC for each procedure to keep low CV value and ensure good reproducibility.

For EH0201 Human IL-6 ELISA Kit, the intra assay CV value is <8%, inter assay CV value is <10%.

Sample ID	Intra assay precision			Inter assay precision		
	1#	2#	3#	1#	2#	3#
n	20	20	20	20	20	20
Mean(pg/ml)	210	113	61	203	105	57
Standard deviation	3.2	1.4	0.6	6.2	3.0	1.3
CV(%)	2.8	3.5	4.1	3.0	4.1	5.6




**Stability:**

The stability of ELISA Kit is determined by the loss rate of activity. The loss rate of this kit is less than 10% within the expiration date under appropriate storage condition.

Standard(n=5)	37°C for 1 month	4°C for 6 months
Average (%)	80	95-100


## Antibodies

FineTest® have 10,000+ antibodies with validated images available for applications like ELISA, WB, IHC, FC...etc, including Rabbit polyclonal & Mouse monoclonal primary antibodies and various tagged secondary antibodies such as HRP, FITC, Biotin...etc.

 <p><b>anti- 6*His, His-Tag antibody</b> Cat NO.:FNab00008 Clone:Monoclonal Application:ELISA, WB, IP, IF fusion protein with his-tag</p>	 <p><b>anti- GFP tag antibody</b> Cat NO.: FNab03433 Clone:Polyclonal Application:ELISA, WB, IP fusion protein with eGFP tag</p>	 <p><b>anti- MBP tag antibody</b> Cat NO.:FNab05042 Clone:Polyclonal Application:ELISA, WB, IP fusion protein with MBP-tag</p>
--	---	---

### anti- IRAK1 antibody


Cat NO.:FNab04382  
Species reactivity: Human, Mouse, Rat  
Application: ELISA, FC, IF, IHC, IP, WB



IF of MCP-7 cells with FITC Human breast cancer tissue MCP-7 cell lysate

### anti- Beclin 1 antibody


Cat NO.:FNab00860  
Species reactivity: Human, Mouse, Rat  
Application: ELISA, IHC, WB, IF



human colorectal cancer tissue IF of Hydroxychloroquine induced MDA-MB-453s cells with FITC Hydroxychloroquine induced MDA-MB-453s cell lysate

## Recombinant Protein

FineTest® protein expression platform has offered nearly over 4500 recombinant proteins expressed by HEK293 and E.Coli for WB, Elisa application. Besides His tag and GST tag, we developed special IF2D1 tag with high-efficient solubility which can greatly improve the protein expression.




P0142 Recombinant Human PAI-1 in SDS Page

P0171 Recombinant Rat IL-1b in SDS Page and ELISA


## IF2D1 Tag

FineTest develops pET28a-IF2D1 expression vector with the powerful fusion tag, which perfectly replaces GST, MBP etc. Furthermore, the molecular weight is lower and the influence on the expression of target proteins is less. Thrombin functions as the enzyme cutting site to facilitate the subsequent cutting.



## HEK293 Cell Expression System

Besides common E.Coli epression system, FineTest also developed HEK293 cell eukaryotic epression system, which can make expressed proteins more naturally active. Meanwhile, with rich epression eperience,our R&D team can customize various protein sequences according to eperimental requirements of customers.



P2750 Recombinant CD147 SDS-PAGE

P2755 Recombinant Human ACE2 SDS-PAGE

Conc. (pg/ml)	OD(450nm)	
	R&D	FineTest®
7.8125	0.079	0.090
15.625	0.308	0.379
31.25	0.532	0.624
62.5	0.858	0.960
125	1.377	1.555
250	2.071	2.297
500	2.853	2.961
1000	3.383	3.322

We conduct the Elisa validation comparison between P2750 Fine Test recombinant CD147 and 972-EMN R&D recombinant -CD147. The results show Fine Test eukaryotic proteins achieve the leading industrial level.


## Service Process



Product Code	Product Name
P2351	Recombinant Mouse Timp2
P2407	Recombinant Human TIMP1
P2536	Recombinant Rabbit TIMP1
P2537	Recombinant Human Timp-2
P2707	Recombinant Human MMP-9
P2750	Recombinant Human EMMPRIN
P2751	Recombinant Human CRP
P3006	Recombinant spike glycoprotein RBD [2019-Novel Coronavirus]
P3018	Recombinant Spike Glycoprotein [2019-Novel Coronavirus]
P2955	Recombinant Human ACE2
...	...

More other items will be available soon.

# Tumor



Tumor-associated antigens (TAA) have elevated levels on tumor cells, but are also expressed at lower levels on healthy cells. Tumor markers are useful in helping to detect, diagnose and manage some types of cancers. The level of present tumor markers can sometimes reflect the stage (extent) of the disease or the patient's prognosis (the likely outcome or the course of disease) for some types of cancers. Many tumor-associated antigens have been identified and targeting them is considered as an important approach in cancer immunotherapy.


FineTest® offers customers a variety of Tumor Research products with high specificity, sensitivity and purity, namely Elisa kits, Antibodies and Proteins. With the high product quality, professional tech support and excellent services, FineTest products sold worldwide are trusted by researchers. Cited academic publications rapidly increased year by year. FineTest always keeps improving ourselves to provide more high-quality products for research use.

## Tumor




### Selected Tumor ELISA Kits

Product Code	Product Name	Range	Product Code	Product Name	Range
EH1256	Human MYOG(Myogenin) ELISA Kit	0.625-40ng/ml	EM0492	Mouse Beclin1(Beclin-1) ELISA Kit	0.156-10ng/ml
EH1278	Human STX1A(Syntaxin-1A) ELISA Kit	0.156-10ng/ml	EM1084	Mouse GPC3(Glypican 3) ELISA Kit	0.156-10ng/ml
EH1484	Human SEMG1(Semenogelin-1) ELISA Kit	15.625-1000pg/ml	EM1196	Mouse LTF/LF(Lactoferrin) ELISA Kit	62.5-4000pg/ml
EH1674	Human OCLN(Occludin) ELISA Kit	0.156-10ng/ml	EM1197	Mouse LZM(Lysozyme) ELISA Kit	31.25-2000ng/ml
EH1915	Human FSCN1(Fascin) ELISA Kit	0.781-50ng/ml	EM1396	Mouse TF(Tissue factor) ELISA Kit	6.25-400pg/ml
EH2944	Human Des(Desmin) ELISA Kit	0.156-10ng/ml	EM1407	Mouse TG(Thyroglobulin) ELISA Kit	125-8000pg/ml
EH2965	Human DSP(Desmoplakin) ELISA Kit	0.313-20ng/ml	ER0593	Rat S100A9(Protein S100-A9) ELISA Kit	0.156-10ng/ml
EH3110	Human GAL8(Galectin 8) ELISA Kit	0.313-20ng/ml	ER0742	Rat AQP-2(Aquaporin 2) ELISA Kit	1.563-100ng/ml
EH3397	Human MUC3(Mucin 3) ELISA Kit	31.25-2000pg/ml	ER0891	Rat Cys-C(Cystatin C) ELISA Kit	15.625-1000pg/ml
EH3649	Human Preptin(Preptin) ELISA Kit	3.125-200pg/ml	ER1233	Rat PCIII(Procollagen III) ELISA Kit	0.156-10ng/ml
EH4017	Human β-actin(Beta actin) ELISA Kit	0.313-20ng/ml	ER1249	Rat PDPN(Podoplanin) ELISA Kit	15.625-1000pg/ml
EH4018	Human ACP(Acid phosphatase) ELISA Kit	0.313-20ng/ml	ER1301	Rat PTX3(Pentraxin 3) ELISA Kit	0.156-10ng/ml



Yinan Zhong et al. Adv. Funct. Mater. 2020, 2007544 tested RhoA in MCF-7 tumors using EH2484



Jaime-Sanchez P, et al. J Immunother Cancer 2020;8:e000528. doi:10.1136/jitc-2020-000528 tested HMGB1 in cell supernatants using EM0382

### Selected Tumor Antibodies

Product Code	Product Name	Reactivity	Application	Product Code	Product Name	Reactivity	Application
FNab00006	anti-4EBP1 antibody	Human	ELISA, WB, IF	FNab00200	anti-AFM antibody	Human	ELISA, WB, IP
FNab00199	anti-AFM antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00473	anti-APC antibody	Human, Mouse	ELISA, WB, IHC
FNab00102	anti-ACPP antibody	Human, Rat	ELISA, WB, IHC	FNab00499	anti-APOD antibody	Human	ELISA, WB, IHC
FNab00183	anti-ADRP antibody	Human	ELISA, WB, IHC	FNab00500	anti-APOD antibody	Human	ELISA, WB
FNab00193	anti-AFAP1L1 antibody	Human, Monkey	ELISA, WB, IHC	FNab00229	anti-AH1 antibody	Human, Mouse	ELISA, WB
FNab00194	anti-AFAP1L2 antibody	Human, Mouse, Rat	ELISA, WB	FNab00278	anti-Albumin antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00196	anti-AFF4 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00753	anti-AXIN2 antibody	Human, Mouse, Rat	ELISA, WB
FNab00197	anti-AFF4 antibody	Human, Mouse, Rat	ELISA, WB, IP	FNab00790	anti-BAGE4 antibody	Human	ELISA, IHC
FNab00522	anti-ARAF antibody	Human, Mouse	ELISA, WB, IHC	FNab00796	anti-BAK antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00761	anti-B2M antibody	Human	ELISA, WB, IHC	FNab00797	anti-BAMBI antibody	Human, Mouse, Rat	ELISA, WB, IP




Ding et al. Cancer Cell Int (2020) 20:125 tested KIF15 in 5-mm sections from formalin-fixedparaffin embedded tissue using FNab04551.

### Selected Tumor Proteins

Product Code	Product Name	Species	Source	Tags	Product Code	Product Name	Species	Source	Tags
P0042	Recombinant Human Insulin	Human	E.Coli	His tag	P0284	Recombinant Human IRF7	Human	E.Coli	His tag
P0119	Recombinant Human trx-1	Human	E.Coli	His tag	P0335	Recombinant Human P105	Human	E.Coli	IF2DI tag
P0142	Recombinant Human PAI-1	Human	E.Coli	His tag	P0336	Recombinant Human P65	Human	E.Coli	IF2DI tag
P0156	Recombinant Human TNFA	Human	E.Coli	His tag	P0337	Recombinant Human P100	Human	E.Coli	IF2DI tag
P0183	Recombinant Human GH	Human	E.Coli	His tag	P0348	Recombinant Human HSPA1A	Human	E.Coli	IF2DI tag
P0197	Recombinant Human Galectin	Human	E.Coli	His tag	P0395	Recombinant Human IDH1	Human	E.Coli	IF2DI tag
P0210	Recombinant Human FABP3	Human	E.Coli	His tag	P0454	Recombinant Human SRC	Human	E.Coli	IF2DI tag
P0227	Recombinant Human RASK	Human	E.Coli	GST tag	P0469	Recombinant Human GFAP	Human	E.Coli	His tag

More FineTest products about Tumor Research are available on our website (www.fn-test.com).

# Metabolism



Metabolism refers to chemical reactions where cells and organisms in their living conditions are catalyzed by enzymes. It is the chemical reaction which happens in cells, including energy metabolism and catabolism in the mitochondrion, biosynthesis of fatty acid and proteins etc. Metabolism is the continuous self-renewal process of the living organism. If the cell metabolism is inhibited, the health of human body will be threatened. Common metabolic diseases include diabetes, heart disease, Alzheimer's disease, obesity, hypertension and arthritis etc. All problems of the disease are attributed to reasons of cells finally. Cell metabolism is the beneficial path to study all kinds of diseases.


FineTest® offers customers a variety of Metabolic Research products with high specificity, sensitivity and purity, namely Elisa kits, Antibodies and Proteins. With the high product quality, professional tech support and excellent services, FineTest products sold worldwide are trusted by researchers. Cited academic publications rapidly increased year by year. FineTest always keeps improving ourselves to provide more high-quality products for research use.

## Metabolism



### Selected Metabolism ELISA Kits

Product Code	Product Name	Range	Product Code	Product Name	Range
EH1034	Human MYH2(Myosin-2) ELISA Kit	78.125-5000pg/ml	EM0395	Mouse Vdr(Vitamin D3 receptor) ELISA Kit	0.313-20ng/ml
EH1146	Human BECN1(Beclin-1) ELISA Kit	0.156-10ng/ml	EM0466	Mouse ACTH(Corticotropin) ELISA Kit	15.625-1000pg/ml
EH0018	Human MPO(Myeloperoxidase) ELISA Kit	1.563-100ng/ml	EM0596	Mouse Csn2(Beta-casein) ELISA Kit	0.781-50ng/ml
EH0023	Human TSP-2(Thrombospondin-2) ELISA Kit	0.313-20ng/ml	EM0609	Mouse Ucn2(Urocortin-2) ELISA Kit	78.125-5000pg/ml
EH0028	Human Activin A ELISA Kit	62.5-4000pg/ml	EM0749	Mouse Smox(Spermine oxidase) ELISA Kit	0.313-20ng/ml
EH0033	Human AR(Androgen receptor) ELISA Kit	0.313-20ng/ml	EM0800	Mouse Nox1(NADPH oxidase 1) ELISA Kit	15.625-1000pg/ml
EH0034	Human ANG(Angiogenin) ELISA Kit	78.125-5000pg/ml	ER0307	Rat Vdr(Vitamin D3 receptor) ELISA Kit	0.313-20ng/ml
EH0049	Human CA9(Carbonic anhydrase 9) ELISA Kit	15.625-1000pg/ml	ER0333	Rat Gzmb(Granzyme B) ELISA Kit	15.625-1000pg/ml
EH0051	Human CTSB(Cathepsin B) ELISA Kit	0.156-10ng/ml	ER0592	Rat S100A8(Protein S100-A8) ELISA Kit	3.125-200ng/ml
EH0052	Human CTSD(Cathepsin D) ELISA Kit	46.875-3000pg/ml	ER1200	Rat NSE(Neuron Specific Enolase) ELISA Kit	78.125-5000pg/ml
EH0054	Human CTSS(Cathepsin S) ELISA Kit	78.125-5000pg/ml	ER1236	Rat PCX(Podocalyxin) ELISA Kit	1.563-100ng/ml
EH0072	Human CD13(Aminopeptidase N) ELISA Kit	31.25-2000pg/ml	ER1408	Rat TRY(Trypsin) ELISA Kit	0.156-10ng/ml



Tristan-Noguero et al. Movement Disorders, 2020 tested APOH, APOD and COL6a3 in CSF samples using EH0531, EH2135 and EH2871.

### Selected Metabolism Antibodies

Product Code	Product Name	Reactivity	Application	Product Code	Product Name	Reactivity	Application
FNab00038	anti- ABCD2 antibody	Human, Mouse, Rat	ELISA, WB, IP	FNab00023	anti- AARSD1 antibody	Human, Mouse, Rat	ELISA, WB
FNab00006	anti- 4EBP1 antibody	Human	ELISA, WB, IF	FNab00025	anti- AATF antibody	Human, Mouse, Rat	ELISA, WB, IF
FNab00065	anti- ACADM antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00039	anti- ABCE1 antibody	Human, Mouse	ELISA, WB
FNab00014	anti- A4GALT antibody	Human, Mouse	ELISA, WB	FNab00027	anti- ABCA2 antibody	Human, Mouse, Rat	ELISA, WB, IF
FNab00015	anti- A4GNT antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00028	anti- ABCA7 antibody	Human, Mouse	ELISA, WB
FNab00016	anti- AACCS antibody	Human, Mouse, Rat	ELISA, WB	FNab00029	anti- ABCA8 antibody	Human, Mouse	ELISA, WB
FNab00022	anti- AARS2 antibody	Human	ELISA, WB, IHC	FNab00030	anti- ABCB10 antibody	Human, Mouse, Rat	ELISA, WB, IP
FNab00018	anti- AADAT antibody	Human, Mouse, Rat	ELISA, IHC, WB	FNab00040	anti- ABCF1 antibody	Human, Mouse, Rat	ELISA, WB, IP
FNab00020	anti- AANAT antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00036	anti- ABCD1 antibody	Human	ELISA, WB
FNab00021	anti- AARE antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00037	anti- ABCD1 antibody	Human	ELISA, WB
FNab00041	anti- ABCF2 antibody	Human, Mouse	ELISA, WB, IHC	FNab00051	anti- ABHD6 antibody	Human, Mouse	ELISA, WB
FNab00043	anti- ABCG4 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00056	anti- ABT1 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00044	anti- ABHD1 antibody	Human, Mouse	ELISA, WB	FNab00058	anti- ACAA1 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00045	anti- ABHD12B antibody	Human, Mouse	ELISA, WB, IHC	FNab00059	anti- ACAA2 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00046	anti- ABHD14A antibody	Human, Mouse, Rat	ELISA, WB	FNab00061	anti- ACAD8 antibody	Human, Mouse, Rat	ELISA, WB
FNab00048	anti- ABHD15 antibody	Human, Mouse	ELISA, WB	FNab00062	anti- ACAD9 antibody	Human, Mouse, Rat	ELISA, WB, IHC

### Selected Metabolism Proteins

Product Code	Product Name	Species	Source	Tags	Product Code	Product Name	Species	Source	Tags
P0042	Recombinant Human Insulin	Human	E.Coli	His tag	P0125	Recombinant Human FABP5	Human	E.Coli	His tag
P0067	Recombinant Human β-Actin	Human	E.Coli	His tag	P0127	Recombinant Human FABP7	Human	E.Coli	IF2DI tag
P0108	Recombinant Human HSPA5	Human	E.Coli	His tag	P0128	Recombinant Human FABP8	Human	E.Coli	IF2DI tag
P0116	Recombinant Human GAPDH	Human	E.Coli	His tag	P0135	Recombinant Human Casp3	Human	E.Coli	IF2DI tag
P0117	Recombinant Human YY1	Human	E.Coli	His tag	P0142	Recombinant Human PAI-1	Human	E.Coli	IF2DI tag
P0119	Recombinant Human trx-1	Human	E.Coli	His tag	P0194	Recombinant Human HMGB1	Human	E.Coli	IF2DI tag
P0122	Recombinant Human FABP2	Human	E.Coli	His tag	P0196	Recombinant Human CYPB	Human	E.Coli	IF2DI tag
P0124	Recombinant Human FABP4	Human	E.Coli	His tag	P0200	Recombinant Human ubc13	Human	E.Coli	His tag

More FineTest products about Metabolism Research are available on our website(www.fn-test.com).

# Cardiovascular

## Cardiovascular

The cardiovascular disease(CVD) is involved in heart or vascular disease, which seriously threatens human beings as a common disease. Several chronic non-infectious diseases most harmful to human beings include hypertension, coronary heart disease, cerebral stroke, diabetes, metabolic syndrome and lipid metabolism disorders etc. These diseases still require a lot of research effort to find the treatment solution. Hence, biomarkers produced during the process of CVD play an important role in diagnosis, risk stratification and estimating prognosis of CVD.


These biomarkers include myocardial necrosis biomarker, myocardial ischemia biomarker, inflammatory reaction biomarker, dyslipidemia biomarker and lipid modification biomarker etc. They can reflect the situation of vulnerable patients, covering vulnerable artery, vulnerable blood and vulnerable myocardium.




FineTest® offers customers a variety of Cardiac Research products with high specificity, sensitivity and purity, namely Elisa kits, Antibodies and Proteins. With the high product quality, professional tech support and excellent services, FineTest products sold worldwide are trusted by researchers. Cited academic publications rapidly increased year by year. FineTest always keeps improving ourselves to provide more high-quality products for research use.

### Selected Cardiovascular ELISA Kits

Product Code	Product Name	Range	Product Code	Product Name	Range
EH0004	Human CLU(Clusterin) ELISA Kit	3.125-200ng/ml	EM0066	Mouse DCN(Decorin) ELISA Kit	7.813-500pg/ml
EH2174	Human APLN(Apelin) ELISA Kit	62.5-4000pg/ml	EM0146	Mouse Neupilin-1 ELISA Kit	0.156-10ng/ml
EH2251	Human WNT2(Protein Wnt-2) ELISA Kit	0.156-10ng/ml	EM0352	Mouse Fbln5(Fibulin-5) ELISA Kit	31.25-2000pg/ml
EH2649	Human AP13(Apelin 13) ELISA Kit	125-8000pg/ml	EM0442	Mouse Txn(Thioedoxin) ELISA Kit	62.5-4000pg/ml
EH2880	Human CPP(Copeptin) ELISA Kit	31.25-2000pg/ml	EM0923	Mouse CHE(Chemerin) ELISA Kit	46.875-3000pg/ml
EH3030	Human F7(Coagulation Factor VII) ELISA Kit	31.25-2000pg/ml	EM1239	Mouse NRG-1(Neuregulin 1) ELISA Kit	78.125-5000pg/ml
EH3055	Human FETUB(Fetuin B) ELISA Kit	15.625-1000pg/ml	ER0046	Rat Lum(Lumican) ELISA Kit	0.625-40ng/ml
EH3065	Human FGα(Fibrinogen Alpha) ELISA Kit	0.469-30ng/ml	ER0074	Rat Afp(Alpha-fetoprotein) ELISA Kit	0.313-20ng/ml
EH3215	Human HBβ(Hemoglobin Beta) ELISA Kit	1.563-100ng/ml	ER0235	Rat Pgf(Placenta growth factor) ELISA Kit	7.813-500pg/ml
EH3443	Human NRG-4(Neuregulin 4) ELISA Kit	0.781-50ng/ml	ER0834	Rat CHEM(Chemerin) ELISA Kit	31.25-2000pg/ml
EH3520	Human PCI(Protein C Inhibitor) ELISA Kit	1.25-80ng/ml	ER0867	Rat CT(Calcitonin) ELISA Kit	15.625-1000pg/ml
EH3847	Human TAGLN(Transgelin) ELISA Kit	0.625-40ng/ml	ER1301	Rat PTX3(Pentraxin 3) ELISA Kit	0.156-10ng/ml




Jun-Neng Roan et al. Journal of Heart and Lung Transplantation (2020), doi: <https://doi.org/10.1016/j.healun.2020.05.017> tested ProT in serum using EH1845



Y. Chen et al. Clinical Science (2016) 130, 2105-2120 doi: 10.1042/CS20160085 tested Troponin I in plasma after UTM treatment using ER0870

### Selected Cardiovascular Antibodies

Product Code	Product Name	Reactivity	Application	Product Code	Product Name	Reactivity	Application
FNab00065	anti-ACADM antibody	Human, Mouse, Rat	ELISA, WB	FNab00138	anti-ADAM17 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00014	anti-A4GALT antibody	Human, Mouse	ELISA, WB, IHC	FNab00140	anti-ADAM19 antibody	Human, Mouse	ELISA, WB, IP
FNab00019	anti-AAMP antibody	Human, Mouse, Rat	ELISA, WB	FNab00144	anti-ADAM8 antibody	Human, Mouse	ELISA, WB, IHC
FNab00031	anti-ABCB6 antibody	Human, Mouse, Rat	ELISA, WB	FNab00173	anti-ADM2 antibody	Human	ELISA, IHC
FNab00059	anti-ACAA2 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00178	anti-ADRA1B antibody	Human, Mouse, Rat	ELISA, WB
FNab00080	anti-ACE antibody	Human, Mouse	ELISA, WB	FNab00180	anti-ADRB2 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00121	anti-ACTN2 antibody	Human, Mouse, Rat	ELISA, WB	FNab00183	anti-ADRP antibody	Human	ELISA, WB, IHC
FNab00137	anti-ADAM12 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00212	anti-AGGF1 antibody	Human, Mouse, Rat	ELISA, WB, IHC



Guo et al. Annals of Translational Medicine, Vol 8, No 20 October 2020 tested iNOS and IL-4 for WB in the rat myocardial tissue using FNab04325

### Selected Cardiovascular Proteins

Product Code	Product Name	Species	Source	Tags	Product Code	Product Name	Species	Source	Tags
P0010	Recombinant Human PCT	Human	E.Coli	His tag	P0149	Recombinant Human IFNG	Human	E.Coli	His tag
P0067	Recombinant Human β-Actin	Human	E.Coli	His tag	P0156	Recombinant Human TNFA	Human	E.Coli	His tag
P0116	Recombinant Human GAPDH	Human	E.Coli	His tag	P0160	Recombinant Human CCL2	Human	E.Coli	His tag
P0122	Recombinant Human FABP2	Human	E.Coli	His tag	P0168	Recombinant Human IL-1b	Human	E.Coli	His tag
P0124	Recombinant Human FABP4	Human	E.Coli	His tag	P0174	Recombinant Human IL-6	Human	E.Coli	His tag
P0125	Recombinant Human FABP5	Human	E.Coli	His tag	P0180	Recombinant Human IL-8	Human	E.Coli	His tag
P0132	Recombinant Human CYGB	Human	E.Coli	His tag	P0190	Recombinant Human leptin	Human	E.Coli	His tag
P0142	Recombinant Human PAI-1	Human	E.Coli	His tag	P0195	Recombinant Human Cofilin	Human	E.Coli	His tag

More FineTest products about Cardiovascular Research are available on our website([www.fn-test.com](http://www.fn-test.com)).

# Neuroscience

Neuroscience is the scientific study of nervous systems. Neuroscience can involve research from many branches of science including neurology, brain science, neurobiology, psychology, computer science, prosthetics, medicine, physics, mathematics, pharmacology, biology etc. Studying the structure and function of nervous system is important to explain how the brain works and diagnosis, prevention or treatment of nervous system disease.

Materials of neurobiological research cover lower animals like drosophila and higher organisms like human, rat, mouse etc. Main methodologies include analysis of nucleic acids and proteins, PCR, immunohistochemistry(IHC), western blot(WB) etc. Common neurobiological markers are neuron axon marker Tau, neuron dendrite marker Drebrin/MAP/SAP102, early neuron marker Tubulin and astrocyte marker GFAP/S100 etc.


FineTest® offers customers a variety of Neurobiological Research products with high specificity, sensitivity and purity, namely Elisa kits, Antibodies and Proteins. With the high product quality, professional tech support and excellent services, FineTest products sold worldwide are trusted by researchers. Cited academic publications rapidly increased year by year. FineTest always keeps improving ourselves to provide more high-quality products for research use.

## Neuroscience




### Selected Neuroscience ELISA Kits

Product Code	Product Name	Range	Product Code	Product Name	Range
EH0050	Human CDH2(Neural Cadherin) ELISA Kit	0.156-10ng/ml	EM0146	Mouse Neuropilin-1 ELISA Kit	0.156-10ng/ml
EH0115	Human DLL1(Delta-like protein 1) ELISA Kit	31.25-2000pg/ml	EM0381	Mouse Tac1(Protachykinin-1) ELISA Kit	31.25-2000pg/ml
EH0128	Human FETUA(Fetuin A) ELISA Kit	0.781-50ng/ml	EM0422	Mouse ApoA1(Apolipoprotein A-I) ELISA Kit	7.813-500ng/ml
EH0145	Human GAL3(Galectin-3) ELISA Kit	0.156-10ng/ml	EM0466	Mouse ACTH(Corticotropin) ELISA Kit	15.625-1000pg/ml
EH0146	Human GAL4(Galectin-4) ELISA Kit	0.313-20ng/ml	EM0565	Mouse Slc30a3(Zinc transporter 3) ELISA Kit	78.125-5000pg/ml
EH0147	Human GAL7(Galectin-7) ELISA Kit	125-8000pg/ml	EM0759	Mouse CD68(Macrosialin) ELISA Kit	0.781-50ng/ml
EH0176	Human IL-12(Interleukin 12) ELISA Kit	15.625-1000pg/ml	ER0090	Rat Fth1(Ferritin heavy chain) ELISA Kit	6.25-400ng/ml
EH0184	Human IL-1α(Interleukin 1 Alpha) ELISA Kit	3.906-250pg/ml	ER0271	Rat Penk(Proenkephalin-A) ELISA Kit	0.156-10ng/ml
EH0185	Human IL-1β(Interleukin 1 Beta) ELISA Kit	3.906-250pg/ml	ER0754	Rat Aβ40(Amyloid Beta 40) ELISA Kit	78.125-5000pg/ml
EH0217	Human Leptin R(Leptin receptor) ELISA Kit	15.625-1000pg/ml	ER0846	Rat CNTN4(Contactin 4) ELISA Kit	0.313-20ng/ml
EH0220	Human LUM(Lumican) ELISA Kit	0.313-20ng/ml	ER1122	Rat LEPR(Leptin Receptor) ELISA Kit	31.25-2000pg/ml
EH0229	Human MK(Midkine) ELISA Kit	15.625-1000pg/ml	ER1375	Rat TFR(Transferrin Receptor) ELISA Kit	0.313-20ng/ml





Daniel Asby et al. Brain, awab094, <https://doi.org/10.1093/brain/awab094> tested Fibrinogen in brain tissue homogenates using EH3057



Ahmed S. Kamel et al. Mol Neurobiol 55, 8188-8202 (2018). <https://doi.org/10.1007/s12035-018-0966-3> tested GLUR2 in hippocampal tissue using ER1901

### Selected Neuroscience Antibodies

Product Code	Product Name	Reactivity	Application	Product Code	Product Name	Reactivity	Application
FNab00027	anti-ABCA2 antibody	Human, Mouse, Rat	ELISA, WB, IF	FNab00178	anti-ADRA1B antibody	Human, Mouse, Rat	ELISA, WB
FNab00028	anti-ABCA7 antibody	Human, Mouse	ELISA, WB	FNab00180	anti-ADRB2 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00033	anti-ABCB7 antibody	Human, Mouse, Rat	ELISA, WB, IP	FNab00229	anti-AHI1 antibody	Human, Mouse	ELISA, WB
FNab00043	anti-ABCG4 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00239	anti-AIP/ARA9 antibody	Human, Mouse	ELISA, WB, IHC
FNab00121	anti-ACTN2 antibody	Human, Mouse, Rat	ELISA, WB	FNab00240	anti-AIPL1 antibody	Human, Mouse, Rat	ELISA, IHC, WB
FNab00138	anti-ADAM17 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00282	anti-ALDH1A1 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00144	anti-ADAM8 antibody	Human, Mouse	ELISA, WB, IHC	FNab00283	anti-ALDH1A1 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00151	anti-ADARB1 antibody	Human, Mouse, Rat	ELISA, WB	FNab00289	anti-ALDH1L1 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00162	anti-ADH1B antibody	Human, Mouse, Rat	ELISA, WB	FNab00309	anti-ALK,CD246 antibody	Human, Mouse, Rat	ELISA, IHC, WB
FNab00174	anti-ADNP antibody	Human, Mouse, Rat	ELISA, WB, IF	FNab00310	anti-ALK,CD246 antibody	Human	ELISA, IHC, WB



Peng et al. Journal of Neuroinflammation (2020) 17:81 tested NLRX1 for immunoprecipitation and WB using FNab05759.

### Selected Neuroscience Proteins

Product Code	Product Name	Species	Source	Tags	Product Code	Product Name	Species	Source	Tags
P0127	Recombinant Human FABP7	Human	E.Coli	His tag	P0197	Recombinant Human Galectin	Human	E.Coli	His tag
P0156	Recombinant Human TNFA	Human	E.Coli	His tag	P0203	Recombinant Human AKT1	Human	E.Coli	GST tag
P0160	Recombinant Human CCL2	Human	E.Coli	His tag	P0221	Recombinant Human GSK3B	Human	E.Coli	His tag
P0168	Recombinant Human IL-1b	Human	E.Coli	His tag	P0225	Recombinant Human JNK	Human	E.Coli	His tag
P0174	Recombinant Human IL-6	Human	E.Coli	His tag	P0251	Recombinant Human LUBAC	Human	E.Coli	His tag
P0183	Recombinant Human GH	Human	E.Coli	His tag	P0255	Recombinant Human PTP1B	Human	E.Coli	His tag
P0194	Recombinant Human HMGB1	Human	E.Coli	His tag	P0262	Recombinant Human FOXO1	Human	E.Coli	His tag
P0195	Recombinant Human Cofilin	Human	E.Coli	His tag	P0267	Recombinant Human ELKS	Human	E.Coli	GST tag

More FineTest products about Neuroscience Research are available on our website([www.fn-test.com](http://www.fn-test.com)).


# Signal Transduction

## Signal Transduction




### Selected Signal Transduction ELISA Kits

Product Code	Product Name	Range	Product Code	Product Name	Range
EH0269	Human RETN(Resistin) ELISA Kit	31.25-2000pg/ml	EM0328	Mouse Prok1(Prokineticin-1) ELISA Kit	6.25-400pg/ml
EH0280	Human TFF1(Trefoil factor 1) ELISA Kit	15.625-1000pg/ml	EM0453	Mouse Hcrt(Orexin) ELISA Kit	7.813-500pg/ml
EH0382	Human FE(Ferritin) ELISA Kit	0.156-10ng/ml	EM0466	Mouse ACTH(Corticotropin) ELISA Kit	15.625-1000pg/ml
EH0388	Human PI(Proinsulin) ELISA Kit	31.25-2000pg/ml	EM0509	Mouse Pyy(Peptide YY) ELISA Kit	15.625-1000pg/ml
EH0419	Human PR-3(Proteinase-3) ELISA Kit	0.156-10ng/ml	EM0846	Mouse APO-E(Apolipoprotein E) ELISA Kit	1.563-100ng/ml
EH0533	Human FST(Follistatin) ELISA Kit	0.25-16ng/ml	EM0927	Mouse C-jun(C-jun) ELISA Kit	0.391-25ng/ml
EH1057	Human HP(Haptoglobin) ELISA Kit	1.563-100ng/ml	ER0371	Rat Agt(Angiotensinogen) ELISA Kit	1.563-100ng/ml
EH2031	Human HPX(Hemopexin) ELISA Kit	1.563-100ng/ml	ER0449	Rat Arg1(Arginase-1) ELISA Kit	2.5-160ng/ml
EH2052	Human KHK(Ketohexokinase) ELISA Kit	0.313-20ng/ml	ER0450	Rat Prdx2(Peroxiredoxin-2) ELISA Kit	0.781-50ng/ml
EH2184	Human KL(Klotho) ELISA Kit	7.813-500pg/ml	ER0602	Rat TRY4(Trypsin-4) ELISA Kit	0.156-10ng/ml
EH3221	Human Hcp3(Hepcidin) ELISA Kit	31.25-2000pg/ml	ER0808	Rat CCK-8(Cholecystokinin 8) ELISA Kit	15.625-1000pg/ml
EH3566	Human PHB(Prohibitin) ELISA Kit	0.313-20ng/ml	ER0895	Rat DAO(Diamine Oxidase) ELISA Kit	3.125-200ng/ml



Vincenzo De Iuliis et al. *Oncotarget*, 2018, Vol. 9, (No. 54), pp: 30289-30303 tested EGF and VEGF in cell supernatants using EH0009 and EH0327




C. Wang, et al. *Phytomedicine* 78 (2020) 153305 tested PGE2 in cell medium using ER1800

Signal transduction refers to the basic biological approach where cells respond to the external stimuli, through membrane or intracellular receptors simulated by extracellular informational molecules, influencing the biological function by the switch of cellular signal transduction system. In the signal transduction process, external signal is recognized, amplified and integrated by cells and then raises the response, e.g. changes in enzyme activity, gene expression or ion channel activity. Signal transduction is involved in many processes, e.g. cell growth, proliferation, healing, differentiation, metabolism and signaling. In the pathological situation, the disorder of cell signal transduction pathway in its single or multiple phases can result in the cell metabolism and dysfunction or abnormal growth. Main cell signal transduction pathways include G protein mediated signal transduction, receptor tyrosine protein kinase(RTPK) signal transduction, Non-RTPK signal transduction, receptor guanylate cyclase signal transduction, nuclear receptor signal transduction. Diseases caused by signal transduction disorders include myocardial infarction, coronary heart disease, diabetes, cerebral infarction, pulmonary infarction and sickle cell anemia etc.

### Selected Signal Transduction Antibodies

Product Code	Product Name	Reactivity	Application	Product Code	Product Name	Reactivity	Application
FNab00069	anti-ACAP1 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00140	anti-ADAM19 antibody	Human, Mouse	ELISA, WB, IP
FNab00070	anti-ACAP2 antibody	Human, Mouse, Rat	ELISA, WB	FNab00141	anti-ADAM2 antibody	Human, Mouse, Rat	ELISA, WB
FNab00071	anti-ACAP3 antibody	Human, Mouse, Rat	ELISA, WB	FNab00142	anti-ADAM28 antibody	Human, Mouse, Rat	ELISA, WB
FNab00086	anti-ACK1 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00144	anti-ADAM8 antibody	Human, Mouse	ELISA, WB, IHC
FNab00089	anti-ACOT11 antibody	Human, Mouse, Rat	ELISA, IHC	FNab00145	anti-ADAMDEC1 antibody	Human	ELISA, WB, IHC
FNab00098	anti-ACP1 antibody	Human, Mouse	ELISA, WB, IHC	FNab00149	anti-ADAP2 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00129	anti-ACVR1B antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00155	anti-ADCY1 antibody	Human	ELISA, WB
FNab00130	anti-ACVR1C antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00157	anti-ADCY6 antibody	Human, Mouse, Rat	ELISA, WB, IP
FNab00137	anti-ADAM12 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00158	anti-ADCY8 antibody	Human, Mouse	ELISA, WB
FNab00138	anti-ADAM17 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00159	anti-Adenosine A1	Human, Mouse	ELISA, WB, IHC



Li C et al. *International Journal of Biochemistry and Cell Biology* (2010), <https://doi.org/10.1016/j.biocel.2018.03.020> tested FADD in brain tissue using FNab02941

### Selected Signal Transduction Proteins

Product Code	Product Name	Species	Source	Tags	Product Code	Product Name	Species	Source	Tags
P0010	Recombinant Human PCT	Human	E.Coli	His tag	P0137	Recombinant Human Uter	Human	E.Coli	His tag
P0042	Recombinant Human Insulin	Human	E.Coli	His tag	P0149	Recombinant Human IFNG	Human	E.Coli	His tag
P0109	Recombinant Human CXCL1	Human	E.Coli	His tag	P0156	Recombinant Human TNFA	Human	E.Coli	His tag
P0122	Recombinant Human FABP2	Human	E.Coli	His tag	P0160	Recombinant Human CCL2	Human	E.Coli	His tag
P0124	Recombinant Human FABP4	Human	E.Coli	His tag	P0164	Recombinant Human CXCL9	Human	E.Coli	His tag
P0125	Recombinant Human FABP5	Human	E.Coli	His tag	P0165	Recombinant Human CCL11	Human	E.Coli	His tag
P0127	Recombinant Human FABP7	Human	E.Coli	His tag	P0168	Recombinant Human IL-1b	Human	E.Coli	His tag
P0128	Recombinant Human FABP8	Human	E.Coli	His tag	P0174	Recombinant Human IL-6	Human	E.Coli	His tag

More FineTest products about Signal Transduction Research are available on our website([www.fn-test.com](http://www.fn-test.com)).

FineTest® offers customers a variety of Signal Transduction Research products with high specificity, sensitivity and purity, namely Elisa kits, Antibodies and Proteins. With the high product quality, professional tech support and excellent services, FineTest products sold worldwide are trusted by researchers. Cited academic publications rapidly increased year by year. FineTest always keeps improving ourselves to provide more high-quality products for research use.




# Stem Cells

## Stem Cells




### Selected Stem Cells ELISA Kits

Product Code	Product Name	Range	Product Code	Product Name	Range
EH0071	Human ENG(Endoglin) ELISA Kit	0.156-10ng/ml	EM0307	Mouse Fst(Follistatin) ELISA Kit	0.313-20ng/ml
EH0116	Human E-Cad(E-Cadherin) ELISA Kit	0.156-10ng/ml	EM0324	Mouse Alcam(CD166 antigen) ELISA Kit	0.156-10ng/ml
EH0124	Human SELE(E-Selectin) ELISA Kit	78.125-5000pg/ml	EM0491	Mouse Cdh3(Cadherin-3) ELISA Kit	0.156-10ng/ml
EH0134	Human FN1(Fibronectin) ELISA Kit	1.563-100ng/ml	EM0639	Mouse Tagln3(Transgelin-3) ELISA Kit	0.313-20ng/ml
EH0162	Human LAP(Leucine Peptidase) ELISA Kit	0.781-50ng/ml	EM0746	Mouse Fabp6(Gastrotropin) ELISA Kit	78.125-5000pg/ml
EH0201	Human IL-6(Interleukin 6) ELISA Kit	4.688-300pg/ml	EM1085	Mouse GPC4(Glypican 4) ELISA Kit	31.25-2000pg/ml
EH0205	Human IL-8(Interleukin 8) ELISA Kit	31.25-2000pg/ml	ER0076	Rat Prl(Prolactin) ELISA Kit	0.781-50ng/ml
EH0216	Human LEP(Leptin) ELISA Kit	31.25-2000pg/ml	ER0162	Rat Sost(Sclerostin) ELISA Kit	31.25-2000pg/ml
EH0248	Human OPN(Osteopontin) ELISA Kit	0.156-10ng/ml	ER0342	Rat Cd44(CD44 antigen) ELISA Kit	0.313-20ng/ml
EH0359	Human AFP(alpha Fetoprotein) ELISA Kit	0.313-20ng/ml	ER0881	Rat CX43(Connexin 43) ELISA Kit	0.156-10ng/ml
EH0388	Human PI(Proinsulin) ELISA Kit	31.25-2000pg/ml	ER1180	Rat NCAD(Neural Cadherin) ELISA Kit	0.156-10ng/ml
EH0422	Human TPO(Thrombopoietin) ELISA Kit	31.25-2000pg/ml	ER1268	Rat PKB(Protein Kinase B) ELISA Kit	0.156-10ng/ml




Wu et al The American Journal of Sports Medicine 1-13 DOI: 10.1177/0363546520907137 tested PRG4 in culture medium using EH1094



Agata Przekora et al. Int. J. Mol. Sci. 2020, 21, 4738; doi:10.3390/ijms21134738 tested bALP in the plasma using EH2691

Stem cells are a class of cells with unlimited or immortal ability to self-renew, capable of producing at least one type of highly differentiated progeny. The two major types present are embryonic stem cells and somatic (adult) stem cells. As seed cells for cell therapy and tissue and organ replacement therapy, stem cells have broad application prospects in the fields of developmental biology, disease model research, drug development and screening, and cell therapy.

Stem cell markers include embryonic stem cell markers like OCT4, SSEA-1; hematopoietic stem cell markers like CD34, CD133; mesenchymal stem cell markers like STRO-1; neural stem cell markers like Nestin, p75NTR.



FineTest® offers customers a variety of Stem Cell Research products with high specificity, sensitivity and purity, namely Elisa kits, Antibodies and Proteins. With the high product quality, professional tech support and excellent services, FineTest products sold worldwide are trusted by researchers. Cited academic publications rapidly increased year by year. FineTest always keeps improving ourselves to provide more high-quality products for research use.

### Selected Stem Cells Antibodies

Product Code	Product Name	Reactivity	Application	Product Code	Product Name	Reactivity	Application
FNab00129	anti- ACVR1B antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00992	anti- BVES antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00130	anti- ACVR1C antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab01065	anti- C1orf31 antibody	Human, Mouse	ELISA, WB, IP
FNab00213	anti- Aggrecan antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab01161	anti- CA2 antibody	Human, Mouse, Rat	ELISA, WB
FNab00278	anti- Albumin antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab01329	anti- CBX2 antibody	Human, Mouse, Rat	ELISA, WB
FNab00280	anti- ALCAM antibody	Human	ELISA, IHC, WB	FNab01416	anti- CD13 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00282	anti- ALDH1A1 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab01417	anti- CD13 antibody	Human	ELISA, WB, IHC
FNab00283	anti- ALDH1A1 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab01419	anti- CD133 antibody	Human, Mouse, Rat	ELISA, WB, FC
FNab00349	anti- ALR antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab01427	anti- CD14 antibody	Human, Mouse	ELISA, WB
FNab00651	anti- ASTL antibody	Human	ELISA, WB	FNab01431	anti- CD151 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00771	anti- B4GALNT1 antibody	Human	ELISA, WB	FNab01438	anti- CD2 antibody	Human, Mouse, Rat	ELISA, IHC, WB
FNab00772	anti- B4GALNT2 antibody	Human	ELISA, WB, IHC	FNab01439	anti- CD2 antibody	Human	ELISA, WB
FNab00917	anti- BMP2 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab01444	anti- CD22 antibody	Human	ELISA, WB, IHC
FNab00919	anti- BMP4 antibody	Human, Mouse, Rat	ELISA, IHC, WB	FNab01449	anti- CD24 antibody	Human, Mouse, Rat	ELISA, WB
FNab00923	anti- BMPR2 antibody	Human, Mouse, Rat	ELISA, WB, IP	FNab01450	anti- CD24 antibody	Human, Mouse	ELISA, WB, IHC
FNab00924	anti- BMPR2 antibody	Human, Mouse, Rat	ELISA, WB, IP	FNab01453	anti- CD27 antibody	Human	ELISA, WB, FC
FNab00969	anti- BSA antibody	Bovine	ELISA, WB	FNab01463	anti- CD31 antibody	Human, Mouse, Rat	ELISA, IHC, WB

### Selected Stem Cells Proteins

Product Code	Product Name	Species	Source	Tags	Product Code	Product Name	Species	Source	Tags
P0042	Recombinant Human Insulin	Human	E.Coli	His tag	P0336	Recombinant Human P65	Human	E.Coli	IF2DI tag
P0122	Recombinant Human FABP2	Human	E.Coli	His tag	P0351	Recombinant Human HIF-1A	Human	E.Coli	His tag
P0124	Recombinant Human FABP4	Human	E.Coli	His tag	P0400	Recombinant Human CAV1	Human	E.Coli	IF2DI tag
P0127	Recombinant Human FABP7	Human	E.Coli	His tag	P0453	Recombinant Human GLUT1	Human	E.Coli	IF2DI tag
P0128	Recombinant Human FABP8	Human	E.Coli	His tag	P0459	Recombinant Human SIRT1	Human	E.Coli	His tag
P0190	Recombinant Human leptin	Human	E.Coli	His tag	P0469	Recombinant Human GFAP	Human	E.Coli	His tag
P0203	Recombinant Human AKT1	Human	E.Coli	GST tag	P0476	Recombinant Human ALPL	Human	E.Coli	IF2DI tag
P0210	Recombinant Human FABP3	Human	E.Coli	His tag	P0483	Recombinant Human BMP2	Human	E.Coli	IF2DI tag

More FineTest products about Stem Cells Research are available on our website(www.fn-test.com).



# Immunology

Immunology is a biomedical science that studies the innate non-specific immune system and the adaptive or specific immune system. In response to disease or infection, the innate immune system provides a rapid, general, non-specific immune response, and if it is unsuccessful in destroying the pathogen, the adaptive immune system then provides a highly specific immune response to target the pathogen. The adaptive immune system is also responsible for the development and maintenance of memory cells, specifically T and B cells, which remember previous aggressors and provide quicker specific immune responses against re-offending pathogens.

Immunological disorders bring harmful reactions and results to the organism, e.g. hypersensitivity reactions, autoimmune diseases and tumour etc. Tumor has been the primary reason resulting in global death. When lesions of tumour appear in the body, special proteins, enzymes and hormones which are hardly seen in healthy conditions will appear in blood and urine. These special substances tested in the medical examination are called tumor markers. Usually, in Japanese medical examination, tumor markers are screened by highly accurate blood and serum test.


FineTest® offers customers a variety of Immunological Research products with high specificity, sensitivity and purity, namely Elisa kits, Antibodies and Proteins. With the high product quality, professional tech support and excellent services, FineTest products sold worldwide are trusted by researchers. Cited academic publications rapidly increased year by year. FineTest always keeps improving ourselves to provide more high-quality products for research use.

## Immunology



### Selected Immunology ELISA Kits

Product Code	Product Name	Range	Product Code	Product Name	Range
EH0001	Human ALCAM(CD166 antigen) ELISA Kit	62.5-4000pg/ml	EM0073	Mouse EPO(Erythropoietin) ELISA Kit	46.875-3000pg/ml
EH0018	Human MPO(Myeloperoxidase) ELISA Kit	1.563-100ng/ml	EM0160	Mouse SELP(P-Selectin) ELISA Kit	15.625-1000pg/ml
EH0023	Human TSP-2(Thrombospondin-2) ELISA Kit	0.313-20ng/ml	EM1146	Mouse IFN-α(Interferon α) ELISA Kit	15.625-1000pg/ml
EH0028	Human Activin A ELISA Kit	62.5-4000pg/ml	EM1197	Mouse LZM(Lysozyme) ELISA Kit	31.25-2000ng/ml
EH0034	Human ANG(Angiogenin) ELISA Kit	78.125-5000pg/ml	EM1248	Mouse OC/BGP(Osteocalcin) ELISA Kit	0.391-25ng/ml
EH0050	Human CDH2(Neural Cadherin) ELISA Kit	0.156-10ng/ml	EM1416	Mouse TP53(Tumor Protein 53) ELISA Kit	62.5-4000pg/ml
EH0051	Human CTSB(Cathepsin B) ELISA Kit	0.156-10ng/ml	ER0006	Rat ADP/Acrp30(Adiponectin) ELISA Kit	15.625-1000pg/ml
EH0070	Human MME(Nepriylsin) ELISA Kit	0.156-10ng/ml	ER0011	Rat TLR7(Toll-Like Receptor 7) ELISA Kit	0.313-20ng/ml
EH0071	Human ENG(Endoglin) ELISA Kit	0.156-10ng/ml	ER0074	Rat Afp(Alpha-fetoprotein) ELISA Kit	0.313-20ng/ml
EH0099	Human CRP(C-Reactive Protein) ELISA Kit	31.25-2000pg/ml	ER0202	Rat TPO(Thrombopoietin) ELISA Kit	31.25-2000pg/ml
EH0116	Human E-Cad(E-Cadherin) ELISA Kit	0.156-10ng/ml	ER0304	Rat C-P(C-Peptide) ELISA Kit	78.125-5000pg/ml
EH0135	Human FCN1(Ficolin-1) ELISA Kit	0.156-10ng/ml	ER1200	Rat NSE(Neuron Specific Enolase) ELISA Kit	78.125-5000pg/ml



Paulina Żelechowska, et al., Cellular Immunology, <https://doi.org/10.1016/j.cellimm.2020.104079> tested IFN-α and IFN-β in cell supernatants using ER1070 and ER1071

### Selected Immunology Antibodies

Product Code	Product Name	Reactivity	Application	Product Code	Product Name	Reactivity	Application
FNab00019	anti-AAMP antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00359	anti-AMH antibody	Human, Mouse	ELISA, WB, IHC
FNab00080	anti-ACE antibody	Human, Mouse	ELISA, WB	FNab00360	anti-AMIGO2 antibody	Human, Mouse, Rat	ELISA, WB
FNab00102	anti-ACPP antibody	Human, rat	ELISA, WB, IHC	FNab00361	anti-AMIGO3 antibody	Human, Mouse, Rat	ELISA, WB
FNab00121	anti-ACTN2 antibody	Human, Mouse, Rat	ELISA, WB	FNab00400	anti-ANGPTL5 antibody	Human	ELISA, WB, IHC
FNab00122	anti-ACTN3 antibody	Human, Mouse	ELISA, WB	FNab00451	anti-ANTXR2 antibody	Human, Mouse, Rat	ELISA, IHC, WB
FNab00138	anti-ADAM17 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00481	anti-APH1A antibody	Human, Mouse, Rat	ELISA, WB
FNab00144	anti-ADAM8 antibody	Human, Mouse	ELISA, WB, IHC	FNab00499	anti-APOD antibody	Human	ELISA, WB, IHC
FNab00151	anti-ADARB1 antibody	Human, Mouse, Rat	ELISA, WB	FNab00500	anti-APOD antibody	Human	ELISA, WB
FNab00209	anti-AGER antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00635	anti-ASGR1 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00212	anti-AGGF1 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00640	anti-ASK1 antibody	Human, Mouse, Rat	ELISA, WB
FNab00280	anti-ALCAM antibody	Human	ELISA, IHC, WB	FNab00677	anti-ATG5 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00282	anti-ALDH1A1 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00678	anti-ATG5 antibody	Human, Mouse, Rat	ELISA, WB
FNab00283	anti-ALDH1A1 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00679	anti-ATG7 antibody	Human, Mouse, Rat	ELISA
FNab00309	anti-ALK,CD246 antibody	Human, Mouse, Rat	ELISA, IHC, WB	FNab00761	anti-B2M antibody	Human	ELISA, WB, IHC
FNab00310	anti-ALK,CD246 antibody	Human	ELISA, IHC, WB	FNab00796	anti-BAK antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00319	anti-Alpha actin antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00835	anti-BCL10 antibody	Human, Mouse, Rat	ELISA, WB, IHC

### Selected Immunology Proteins

Product Code	Product Name	Species	Source	Tags	Product Code	Product Name	Species	Source	Tags
P0010	Recombinant Human PCT	Human	E.Coli	His tag	P0168	Recombinant Human IL-1b	Human	E.Coli	His tag
P0042	Recombinant Human Insulin	Human	E.Coli	His tag	P0174	Recombinant Human IL-6	Human	E.Coli	His tag
P0109	Recombinant Human CXCL1	Human	E.Coli	His tag	P0180	Recombinant Human IL-8	Human	E.Coli	His tag
P0149	Recombinant Human IFNG	Human	E.Coli	His tag	P0183	Recombinant Human GH	Human	E.Coli	His tag
P0156	Recombinant Human TNFA	Human	E.Coli	His tag	P0190	Recombinant Human leptin	Human	E.Coli	His tag
P0160	Recombinant Human CCL2	Human	E.Coli	His tag	P0195	Recombinant Human Cofilin	Human	E.Coli	His tag
P0164	Recombinant Human CXCL9	Human	E.Coli	GST tag	P0197	Recombinant Human Galectin	Human	E.Coli	His tag
P0165	Recombinant Human CCL11	Human	E.Coli	His tag	P0200	Recombinant Human ubc13	Human	E.Coli	GST tag

More FineTest products about Immunology Research are available on our website([www.fn-test.com](http://www.fn-test.com)).


# Cell Division and Proliferation

## Cell Division and Proliferation



### Selected Cell Division and Proliferation ELISA Kits

Product Code	Product Name	Range	Product Code	Product Name	Range
EH0706	Human CORO1A(Coronin-1A) ELISA Kit	0.156-10ng/ml	EM0169	Mouse SAA(Serum amyloid A) ELISA Kit	78.125-5000pg/ml
EH1332	Human C19orf80(Betatrophin) ELISA Kit	46.875-3000pg/ml	EM0327	Mouse Coro1a(Coronin-1A) ELISA Kit	0.156-10ng/ml
EH1468	Human HRAS(GTPase HRas) ELISA Kit	0.156-10ng/ml	EM1051	Mouse GAL1(Galectin 1) ELISA Kit	0.781-50ng/ml
EH1617	Human NDRG1(Protein NDRG1) ELISA Kit	0.625-40ng/ml	EM1248	Mouse OC/BGP(Osteocalcin) ELISA Kit	0.391-25ng/ml
EH2013	Human EREG(Proepiregulin) ELISA Kit	15.625-1000pg/ml	EM1285	Mouse PHB(Prohibitin) ELISA Kit	0.391-25ng/ml
EH2186	Human SEPT5(Septin-5) ELISA Kit	0.781-50ng/ml	EM1416	Mouse TP53(Tumor Protein 53) ELISA Kit	62.5-4000pg/ml
EH2487	Human AXIN2(Axin-2) ELISA Kit	0.313-20ng/ml	ER0018	Rat EGF(Epidermal Growth Factor) ELISA Kit	15.625-1000pg/ml
EH2752	Human CALD(Caldesmon) ELISA Kit	0.313-20ng/ml	ER0424	Rat Esr1(Estrogen receptor) ELISA Kit	0.313-20ng/ml
EH2776	Human CCNB(Cyclin-B) ELISA Kit	0.313-20ng/ml	ER0812	Rat CCND3(Cyclin D3) ELISA Kit	31.25-2000pg/ml
EH3072	Human FLOT2(Flotillin 2) ELISA Kit	0.625-40ng/ml	ER0972	Rat GAL3(Galectin 3) ELISA Kit	3.125-200pg/ml
EH3103	Human GAL12(Galectin 12) ELISA Kit	0.313-20ng/ml	ER1274	Rat PKCa(Protein Kinase Ca) ELISA Kit	0.313-20ng/ml
EH3104	Human GAL14(Galectin 14) ELISA Kit	0.313-20ng/ml	ER1439	Rat $\beta$ -actin(Beta actin) ELISA Kit	0.313-20ng/ml




Dai et al., Cell (2019), <https://doi.org/10.1016/j.cell.2019.01.016> tested PGE2 in cell lysates using EH4233

Ni S et al. Cell Prolif. 2020;00:e12882. <https://doi.org/10.1111/cpr.12882> tested NOS2, TRX1 and GSR in cell supernatant using EM0272, EM0442 and EM0720.

Cell proliferation is the process that results in an increase of the number of cells, and is defined by the balance between cell divisions and cell loss through cell death or differentiation. Cell division is the process by which a parent cell divides into two or more daughter cells. Cell division usually occurs as part of a larger cell cycle. It's known that cell proliferation is closely relevant to tumorigenesis and tumor development. Because tumor is a kind of cell cycle disease, the proliferation research is very important.

### Selected Cell Division and Proliferation Antibodies

Product Code	Product Name	Reactivity	Application	Product Code	Product Name	Reactivity	Application
FNab00119	anti- ACTL7A antibody	Human, Mouse, Rat	ELISA, WB, IF	FNab00473	anti- APC antibody	Human, Mouse	ELISA, WB, IHC
FNab00120	anti- ACTL7B antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00475	anti- APC11 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00127	anti- ACTR3B antibody	Human, Mouse, Rat	ELISA, WB	FNab00476	anti- APC2 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00128	anti- ACTRT2 antibody	Human, Mouse, Rat	ELISA, WB	FNab00477	anti- APC4 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00242	anti- Ajuba antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00581	anti- ARL8A antibody	Human, Mouse, Rat	ELISA, WB, IP
FNab00258	anti- AKAP8L antibody	Human	ELISA, WB, IF	FNab00582	anti- ARL8B antibody	Human, Mouse	ELISA, IHC, WB
FNab00348	anti- ALPL2 antibody	Human, Mouse, Rat	ELISA, IHC	FNab00597	anti- ARP2 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00350	anti- ALS2 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00598	anti- ARP3 antibody	Human, Mouse, Rat	ELISA, WB, IP



Taeseung Lee, et al., Molecular and Cellular Endocrinology, <https://doi.org/10.1016/j.mce.2019.03.001> tested Asprosin for Immunoprecipitation from MIN6 cells using FNab09797.

### Selected Cell Division and Proliferation Proteins

Product Code	Product Name	Species	Source	Tags	Product Code	Product Name	Species	Source	Tags
P0067	Recombinant Human $\beta$ -Actin	Human	E.Coli	His tag	P0328	Recombinant Human LKB1	Human	E.Coli	IF2DI tag
P0221	Recombinant Human GSK3B	Human	E.Coli	His tag	P0336	Recombinant Human P65	Human	E.Coli	IF2DI tag
P0227	Recombinant Human RASK	Human	E.Coli	GST tag	P0401	Recombinant Human Irisin	Human	E.Coli	IF2DI tag
P0299	Recombinant Human GDF2	Human	E.Coli	IF2DI tag	P0475	Recombinant Human Actc1	Human	E.Coli	IF2DI tag
P0302	Recombinant Human GDF8	Human	E.Coli	IF2DI tag	P0476	Recombinant Human ALPL	Human	E.Coli	IF2DI tag
P0315	Recombinant Human ERK2	Human	E.Coli	IF2DI tag	P0500	Recombinant Human BRAF	Human	E.Coli	IF2DI tag
P0323	Recombinant Human PCAF	Human	E.Coli	IF2DI tag	P0589	Recombinant Human CDK4	Human	E.Coli	IF2DI tag
P0325	Recombinant Human CK2A1	Human	E.Coli	IF2DI tag	P0616	Recombinant Human MET	Human	E.Coli	IF2DI tag

More FineTest products about Cell Division and Proliferation Research are available on our website([www.fn-test.com](http://www.fn-test.com)).

FineTest® offers customers a variety of Cell Division and Proliferation Research products with high specificity, sensitivity and purity, namely Elisa kits, Antibodies and Proteins. With the high product quality, professional tech support and excellent services, FineTest products sold worldwide are trusted by researchers. Cited academic publications rapidly increased year by year. FineTest always keeps improving ourselves to provide more high-quality products for research use.

# Developmental Biology


## Developmental Biology




### Selected Developmental Biology ELISA Kits

Product Code	Product Name	Range	Product Code	Product Name	Range
EH0119	Human ES(Endostatin) ELISA Kit	0.156-10ng/ml	EM0066	Mouse DCN(Decorin) ELISA Kit	7.813-500pg/ml
EH0128	Human FETUA(Fetuin A) ELISA Kit	0.781-50ng/ml	EM0146	Mouse Neupilin-1 ELISA Kit	0.156-10ng/ml
EH0196	Human IL-3(Interleukin 3) ELISA Kit	15.625-1000pg/ml	EM0290	Mouse POSTN/OSF-2(Periostin) ELISA Kit	0.156-10ng/ml
EH0217	Human Leptin R(Leptin receptor) ELISA Kit	15.625-1000pg/ml	EM0618	Mouse Mgp(Matrix Gla protein) ELISA Kit	78.125-5000pg/ml
EH0225	Human C-ERC(Mesothelin) ELISA Kit	62.5-4000pg/ml	EM0705	Mouse Eda(Ectodysplasin-A) ELISA Kit	0.156-10ng/ml
EH0229	Human MK(Midkine) ELISA Kit	15.625-1000pg/ml	EM0757	Mouse Amot(Angiomotin) ELISA Kit	78.125-5000pg/ml
EH0246	Human NT-4(Neurotrophin-4) ELISA Kit	93.75-6000pg/ml	ER0713	Rat ADM(Adrenomedullin) ELISA Kit	15.625-1000pg/ml
EH0249	Human OSM(Oncostatin-M) ELISA Kit	15.625-1000pg/ml	ER0828	Rat CFH(Complement Factor H) ELISA Kit	15.625-1000ng/ml
EH0255	Human POSTN/OSF2(Periostin) ELISA Kit	0.156-10ng/ml	ER0993	Rat GH(Growth Hormone) ELISA Kit	0.313-20ng/ml
EH0268	Human REN(Renin) ELISA Kit	31.25-2000pg/ml	ER1235	Rat PCT(Procalcitonin) ELISA Kit	15.625-1000pg/ml
EH0279	Human CLEC3B(Tetranectin) ELISA Kit	0.781-50ng/ml	ER1274	Rat PKCa(Protein Kinase Ca) ELISA Kit	0.313-20ng/ml
EH0282	Human TFF3(Trefoil factor 3) ELISA Kit	78.125-5000pg/ml	ER1331	Rat SEMA7A(Semaphorin-7A) ELISA Kit	0.313-20ng/ml

Developmental biology is at the core of all biology. It deals with the process by which the genes in the fertilized egg control cell behavior in the embryo and so determine its pattern, its form, and much of its behavior. Understanding how embryos develop can help to understand developmental abnormalities and other medical challenges, such as cancer. Developmental biology mainly studies the mechanism of multicellular organisms, covering the life process of germ cell from germination, fertilization and embryonic development to its senescence and death.




Wang et al. Cell Death and Disease (2021) 12:212 tested Adiponectin, TNF-α, IL-1β, IL-6 in the serum using EM0001, EM0183, EM0109 and EM0121.



### Selected Developmental Biology Antibodies

Product Code	Product Name	Reactivity	Application	Product Code	Product Name	Reactivity	Application
FNab00019	anti-AAMP antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00178	anti-ADRA1B antibody	Human, Mouse, Rat	ELISA, WB
FNab00034	anti-ABC9 antibody	Human, Mouse	ELISA, WB, IF	FNab00186	anti-Advillin antibody	Human, Mouse, Rat	ELISA, WB
FNab00054	anti-ABLIM1 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00209	anti-AGER antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00055	anti-ABLIM2 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00212	anti-AGGF1 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00138	anti-ADAM17 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00213	anti-Aggregan antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00140	anti-ADAM19 antibody	Human, Mouse	ELISA, WB, IP	FNab00229	anti-AHI1 antibody	Human, Mouse	ELISA, WB
FNab00141	anti-ADAM2 antibody	Human, Mouse, Rat	ELISA, WB	FNab00230	anti-AHNAK antibody	Human	ELISA, IF
FNab00142	anti-ADAM28 antibody	Human, Mouse, Rat	ELISA, WB	FNab00234	anti-AIDA antibody	Human	ELISA, WB, IP



Chen et al. AGING 2020, Vol. 12, No. 10 tested KIF15 for WB, IHC and IF in cell using FNab04551

FineTest® offers customers a variety of Developmental Biology Research products with high specificity, sensitivity and purity, namely Elisa kits, Antibodies and Proteins. With the high product quality, professional tech support and excellent services, FineTest products sold worldwide are trusted by researchers. Cited academic publications rapidly increased year by year. FineTest always keeps improving ourselves to provide more high-quality products for research use.

### Selected Developmental Biology Proteins

Product Code	Product Name	Species	Source	Tags	Product Code	Product Name	Species	Source	Tags
P0010	Recombinant Human PCT	Human	E.Coli	His tag	P0183	Recombinant Human GH	Human	E.Coli	His tag
P0109	Recombinant Human CXCL1	Human	E.Coli	His tag	P0203	Recombinant Human AKT1	Human	E.Coli	GST tag
P0122	Recombinant Human FABP2	Human	E.Coli	His tag	P0207	Recombinant Human CXCL1	Human	E.Coli	His tag
P0124	Recombinant Human FABP4	Human	E.Coli	His tag	P0210	Recombinant Human FABP3	Human	E.Coli	His tag
P0125	Recombinant Human FABP5	Human	E.Coli	His tag	P0221	Recombinant Human GSK3B	Human	E.Coli	His tag
P0127	Recombinant Human FABP7	Human	E.Coli	His tag	P0234	Recombinant Human IRAK1	Human	E.Coli	GST tag
P0128	Recombinant Human FABP8	Human	E.Coli	His tag	P0267	Recombinant Human ELKS	Human	E.Coli	GST tag
P0180	Recombinant Human IL-8	Human	E.Coli	His tag	P0299	Recombinant Human GDF2	Human	E.Coli	IF2DI tag

More FineTest products about Developmental biology Research are available on our website(www.fn-test.com).


# Epigenetics

## Epigenetics




### Selected Epigenetics ELISA Kits

Product Code	Product Name	Range	Product Code	Product Name	Range
EH0099	Human CRP(C-Reactive Protein) ELISA Kit	31.25-2000pg/ml	EM0637	Mouse H2afx(Histone H2A.x) ELISA Kit	31.25-2000pg/ml
EH0115	Human DLL1(Delta-like protein 1) ELISA Kit	31.25-2000pg/ml	EM0861	Mouse ATXN1(Ataxin 1) ELISA Kit	15.625-1000pg/ml
EH0571	Human CTNNβ1(Catenin beta-1) ELISA Kit	0.156-10ng/ml	EM0957	Mouse CTNNβ1(Catenin, Beta 1) ELISA Kit	0.156-10ng/ml
EH0646	Human C4A(Complement C4-A) ELISA Kit	1.563-100ng/ml	EM1117	Mouse histon-H2b(histon-H2b) ELISA Kit	0.391-25ng/ml
EH0661	Human C3(Complement C3) ELISA Kit	7.813-500ng/ml	EM1133	Mouse HSF1(Heat Shock Factor 1) ELISA Kit	0.781-50ng/ml
EH0680	Human C5(Complement C5) ELISA Kit	3.125-200ng/ml	EM1416	Mouse TP53(Tumor Protein 53) ELISA Kit	62.5-4000pg/ml
EH0886	Human MAX(Protein max) ELISA Kit	0.156-10ng/ml	ER0390	Rat Hdac1(Histone deacetylase 1) ELISA Kit	0.156-10ng/ml
EH0115	Human DLL1(Delta-like protein 1) ELISA Kit	31.25-2000pg/ml	ER0452	Rat Dnase1(Deoxyribonuclease-1) ELISA Kit	0.391-25ng/ml
EH0571	Human CTNNβ1(Catenin beta-1) ELISA Kit	0.156-10ng/ml	ER0459	Rat Maf(Transcription factor Maf) ELISA Kit	0.156-10ng/ml
EH0661	Human C3(Complement C3) ELISA Kit	7.813-500ng/ml	ER1033	Rat HDAC2(Histone Deacetylase 2) ELISA Kit	0.156-10ng/ml
EH0886	Human MAX(Protein max) ELISA Kit	0.156-10ng/ml	ER1186	Rat NF-κB(Nuclear Factor Kappa B) ELISA Kit	0.156-10ng/ml
EH1256	Human MYOG(Myogenin) ELISA Kit	0.625-40ng/ml	ER1640	Rat HD3(Histone Deacetylase 3) ELISA Kit	0.156-10ng/ml



KWAK et al. The FASEB Journal. 2020;00:1-16. <https://doi.org/10.1096/fj.202000925R> tested CHI3L1 in the plasma using EH0093




Shu Yih Chew et al. Scientific Reports (2019) 9:2843 <https://doi.org/10.1038/s41598-019-39117-1> tested MPO in kidney homogenates using EM0010.

Epigenetics is a branch of genetics that studies heritable changes in gene expression without changes in nucleotide sequence. There are many epigenetic phenomena, including DNA methylation, genomic imprinting, maternal effects, gene silencing, nucleolar dominance, dormant transposon activation and RNA editing. With the development of cell and molecular biology and numerous reported research achievements on human genome project, people gradually have a clear recognition on the mechanism of tumorigenesis and tumor development. The research on tumor markers has been greatly extended and updated. Especially, the development of genomics and epigenetics greatly promotes the discovery of new tumor markers.

### Selected Epigenetics Antibodies

Product Code	Product Name	Reactivity	Application	Product Code	Product Name	Reactivity	Application
FNab00019	anti-AAMP antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00178	anti-ADRA1B antibody	Human, Mouse, Rat	ELISA, WB
FNab00034	anti-ABC9 antibody	Human, Mouse	ELISA, WB, IF	FNab00186	anti-Advillin antibody	Human, Mouse, Rat	ELISA, WB
FNab00054	anti-ABLIM1 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00209	anti-AGER antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00055	anti-ABLIM2 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00212	anti-AGGF1 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00138	anti-ADAM17 antibody	Human, Mouse, Rat	ELISA, WB, IHC	FNab00213	anti-Aggregan antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00140	anti-ADAM19 antibody	Human, Mouse	ELISA, WB, IP	FNab00229	anti-AHII antibody	Human, Mouse	ELISA, WB
FNab00141	anti-ADAM2 antibody	Human, Mouse, Rat	ELISA, WB	FNab00230	anti-AHNAK antibody	Human	ELISA, IF
FNab00142	anti-ADAM28 antibody	Human, Mouse, Rat	ELISA, WB	FNab00234	anti-AIDA antibody	Human	ELISA, WB, IP



Ramóna Pap et al. Int. J. Mol. Sci. 2020, 21, 3175; doi:10.3390/ijms21093175 tested progesterone receptor, fractalkine receptor, actin receptor 1B, matrix-metalloproteinase 2, matrix-metalloproteinase 9, SRC-1, Sox17, and BMP2 for WB using FNab09774, FNab02088, FNab00129, FNab05238, FNab05247, FNab08216, FNab08124 and FNab00917.

### Selected Epigenetics Proteins

Product Code	Product Name	Species	Source	Tags	Product Code	Product Name	Species	Source	Tags
P0194	Recombinant Human HMGB1	Human	E.Coli	His tag	P0303	Recombinant Human TRAF2	Human	E.Coli	IF2DI tag
P0200	Recombinant Human ubc13	Human	E.Coli	GST tag	P0304	Recombinant Human CYLD	Human	E.Coli	IF2DI tag
P0243	Recombinant Human IKKα	Human	E.Coli	GST tag	P0305	Recombinant Human FYN	Human	E.Coli	IF2DI tag
P0249	Recombinant Human HDAC1	Human	E.Coli	GST tag	P0310	Recombinant Human HDAC2	Human	E.Coli	IF2DI tag
P0262	Recombinant Human FOXO1	Human	E.Coli	His tag	P0316	Recombinant Human HDAC1	Human	E.Coli	IF2DI tag
P0284	Recombinant Human IRF7	Human	E.Coli	His tag	P0323	Recombinant Human PCAF	Human	E.Coli	IF2DI tag

More FineTest products about Epigenetics Research are available on our website([www.fn-test.com](http://www.fn-test.com)).

FineTest® offers customers a variety of Epigenetic Research products with high specificity, sensitivity and purity, namely Elisa kits, Antibodies and Proteins. With the high product quality, professional tech support and excellent services, FineTest products sold worldwide are trusted by researchers. Cited academic publications rapidly increased year by year. FineTest always keeps improving ourselves to provide more high-quality products for research use.

## Hot Selling ELISA Kits

Product Code	Product Name	Range	Product Code	Product Name	Range
EH0011	Human IL-18(Interleukin 18) ELISA Kit	15.625-1000pg/ml	EM0260	Mouse INS(Insulin) ELISA Kit	78.125-5000pg/ml
EH0110	Human CST3(Cystatin C) ELISA Kit	0.313-20ng/ml	EM0466	Mouse ACTH(Corticotropin) ELISA Kit	15.625-1000pg/ml
EH0143	Human FUR(Furin) ELISA Kit	15.625-1000pg/ml	EM0468	Mouse F11(Coagulation factor XI) ELISA Kit	3.125-200ng/ml
EH0164	Human IFN-γ(Interferon gamma) ELISA Kit	15.625-1000pg/ml	EM0950	Mouse Cr(Creatinine) ELISA Kit	1.25-80nmol/ml
EH0173	Human IL-10(Interleukin 10) ELISA Kit	7.813-500pg/ml	EM0979	Mouse D2D(D-Dimer) ELISA Kit	7.813-500ng/ml
EH0174	Human IL-11(Interleukin 11) ELISA Kit	31.25-2000pg/ml	EM1128	Mouse HO-1(Heme Oxygenase 1) ELISA Kit	15.625-1000pg/ml
EH0176	Human IL-12(Interleukin 12) ELISA Kit	15.625-1000pg/ml	EM1466	Mouse CTNI(Cardiac Troponin-I) ELISA Kit	7.813-500pg/ml
EH0185	Human IL-1β(Interleukin 1 Beta) ELISA Kit	3.906-250pg/ml	EM1501	Mouse Estrogen ELISA Kit	15.625-1000pg/ml
EH0199	Human IL-4(Interleukin 4) ELISA Kit	31.25-2000pg/ml	ER0042	Rat IL-6(Interleukin-6) ELISA Kit	62.5-4000pg/ml
EH0201	Human IL-6(Interleukin 6) ELISA Kit	4.688-300pg/ml	ER0371	Rat Agt(Angiotensinogen) ELISA Kit	1.563-100ng/ml
EH0205	Human IL-8(Interleukin 8) ELISA Kit	31.25-2000pg/ml	ER0848	Rat COL1(Collagen Type I) ELISA Kit	0.156-10ng/ml
EH0216	Human LEP(Leptin) ELISA Kit	31.25-2000pg/ml	ER0859	Rat CORT(Corticosterone) ELISA Kit	2.813-180ng/ml
EH0220	Human LUM(Lumican) ELISA Kit	0.313-20ng/ml	ER0895	Rat DAO(Diamine Oxidase) ELISA Kit	3.125-200ng/ml
EH0247	Human OPG(Osteoprotegerin) ELISA Kit	62.5-4000pg/ml	ER0972	Rat GAL3(Galectin 3) ELISA Kit	3.125-200pg/ml
EH0248	Human OPN(Osteopontin) ELISA Kit	0.156-10ng/ml	ER1036	Rat Hepc(Hepcidin) ELISA Kit	15.625-1000pg/ml
EH0378	Human D2D(D-Dimer) ELISA Kit	78.125-5000pg/ml	ER1094	Rat IL-1β(Interleukin 1 Beta) ELISA Kit	31.25-2000pg/ml
EH0397	Human LH(Luteinizing Hormone) ELISA Kit	0.938-60mIU/ml	ER1113	Rat INS(Insulin) ELISA Kit	78.125-5000pg/ml
EH0406	Human MUC-1(Mucin-1) ELISA Kit	0.156-10ng/ml	ER1123	Rat LH(Luteinizing Hormone) ELISA Kit	1.563-100mIU/ml
EH0417	Human IgG(Immunoglobulin G) ELISA Kit	1.563-100ng/ml	ER1432	Rat α1-AGP(α1-Acid Glycoprotein) ELISA Kit	6.25-400ng/ml
EH0468	Human 5-NT(5-Nucleotidase) ELISA Kit	0.313-20ng/ml	ER1462	Rat T(Testosterone) ELISA Kit	0.313-20ng/ml
EH0641	Human COR(Cortisol) ELISA Kit	0.391-25ng/ml	ER1505	Rat Estrogen ELISA Kit	15.625-1000pg/ml
EH0680	Human C5(Complement C5) ELISA Kit	3.125-200ng/ml	EU0253	ST/5-HT(5-hydroxytryptamine) ELISA Kit	1.563-100ng/ml
EH0813	Human CTSK(Cathepsin K) ELISA Kit	0.156-10ng/ml	EU0392	DA(Dopamine) ELISA Kit	1.563-100ng/ml
EH1025	Human TTR(Transferrin) ELISA Kit	15.625-1000ng/ml	EU2539	VB12(Vitamin B12) ELISA Kit	0.781-50ng/ml
EH1042	Human PODXL(Podocalyxin) ELISA Kit	0.156-10ng/ml	EU2554	PGE2(Prostaglandin E2) ELISA Kit	31.25-2000pg/ml
EH1164	Human WNT5A(Protein Wnt-5a) ELISA Kit	0.156-10ng/ml	EU2556	HA(Hyaluronic Acid) ELISA Kit	1.563-100ng/ml
EH1353	Human FLG(Filaggrin) ELISA Kit	0.313-20ng/ml	EU2560	3-NT(3-Nitrotyrosine) ELISA Kit	1.563-100ng/ml
EH1845	Human PTMA(Prothymosin alpha) ELISA Kit	0.156-10ng/ml	EU2563	EPI(Epinephrine/Adrenaline) ELISA Kit	7.813-500pg/ml
EH1927	Human Opi(Opiorphin) ELISA Kit	0.156-10ng/ml	EU2577	MDA(Malondialdehyde) ELISA Kit	7.813-500ng/ml
EH2017	Human PROM1(Prominin-1) ELISA Kit	0.156-10ng/ml	EH4184	Human Estrogen ELISA Kit	15.625-1000pg/ml
EH2080	Human 4-HNE(4-Hydroxynonenal) ELISA Kit	31.25-2000pg/ml	EH4250	Human 3MH(3-methylhistidine) ELISA Kit	6.25-400nmol/ml
EH2119	Human GREM1(Gremelin-1) ELISA Kit	0.156-10ng/ml	EH4320	Human OVA(Ovalbumin) ELISA Kit	0.781-50ng/ml
EH2184	Human KL(Klotho) ELISA Kit	7.813-500pg/ml	EH4704	Human ProCollagen II ELISA Kit	0.156-10ng/ml
EH2381	Human EPX(Eosinophil peroxidase) ELISA Kit	1.563-100ng/ml	EH4810	Human pro-ET1(Pro-endothelin 1) ELISA Kit	15.625-1000pg/ml
EH2552	Human his(Histamine) ELISA Kit	1.563-100ng/ml	EH7533	Human COLI(Collagen Type I) ELISA Kit	0.313-20ng/ml
EH2593	Human ADP(Adiponectin) ELISA Kit	1.563-100ng/ml	EM1592	Mouse IL-8(Interleukin 8) ELISA Kit	15.625-1000pg/ml
EH2618	Human ALP(Alkaline Phosphatase) ELISA Kit	3.125-200ng/ml	EM1600	Mouse FUR(Furin) ELISA Kit	46.875-3000pg/ml
EH2627	Human ANG II(Angiotensin II) ELISA Kit	31.25-2000pg/ml	EM1620	Mouse Calprotectin(Calprotectin) ELISA Kit	7.813-500ng/ml
EH2665	Human ARO(Aromatase) ELISA Kit	0.156-10ng/ml	EM1707	Mouse Asprosin(Asprosin) ELISA Kit	0.313-20ng/ml
EH2685	Human Aβ42(Amyloid Beta 42) ELISA Kit	4.688-300pg/ml	EM1710	Mouse Mcpt4(Mast cell protease 4) ELISA Kit	46.875-3000pg/ml
EH2880	Human CPP(Copeptin) ELISA Kit	31.25-2000pg/ml	ER1070	Rat IFN-α(Interferon α) ELISA Kit	15.625-1000pg/ml
EH3057	Human FG(Fibrinogen) ELISA Kit	1.563-100ng/ml	EU0192	PAF(Platelet Activating Factor) ELISA Kit	0.156-10ng/ml
EH3266	Human IL-13(Interleukin 13) ELISA Kit	15.625-1000pg/ml	EU0199	MT(Melatonin) ELISA Kit	7.813-500pg/ml
EH3267	Human IL-17A(Interleukin 17 A) ELISA Kit	31.25-2000pg/ml	EU0382	HIS(Histamine) ELISA Kit	1.563-100ng/ml
EH3273	Human IL-35(Interleukin 35) ELISA Kit	15.625-1000pg/ml	EU0187	4-HNE(4-Hydroxynonenal) ELISA Kit	31.25-2000pg/ml
EH3301	Human LTA(lipoteichoic acids) ELISA Kit	7.813-500pg/ml	EU0188	KYN(Kynurenine) ELISA Kit	7.813-500pmol/ml
EH3344	Human MT(Melatonin) ELISA Kit	7.813-500pg/ml	EU3118	PD1(Protectin D1) ELISA Kit	0.156-10ng/ml
EH3413	Human Neopterin ELISA Kit	0.156-10ng/ml	ER1878	Rat MDA(Malonaldehyde) ELISA Kit	7.813-500ng/ml
EH3468	Human OC/BGP(Osteocalcin) ELISA Kit	1.25-80ng/ml	EM1905	Mouse Lum(Lumican) ELISA Kit	0.625-40ng/ml
EH3521	Human PCIII(Procollagen III) ELISA Kit	78.125-5000pg/ml	EH4940	Human anti-SARS-CoV2(N) IgM ELISA Kit	0.781-50ng/ml
EH3818	Human sSELP(Soluble P-Selectin) ELISA Kit	15.625-1000pg/ml	EH4397	Human anti-SARS-CoV2(N) IgG ELISA Kit	3.906-250ng/ml
EH3906	Human TPS(Tryptase) ELISA Kit	0.156-10ng/ml	EH4702	Human Irisin ELISA Kit	1.563-100ng/ml
EH4010	Human HS(Heparan Sulfate) ELISA Kit	78.125-5000ng/ml	EH3655	Human pro-IL1β ELISA Kit	31.25-2000pg/ml
EH4011	Human HCY(Homocysteine) ELISA Kit	7.813-500pmol/ml	EU0403	T3(Triiodothyronine) ELISA Kit	0.156-10ng/ml
EH4114	Human Presepsin(Presepsin) ELISA Kit	0.156-10ng/ml	EU0380	PREG(Pregnenolone) ELISA Kit	15.625-1000pg/ml
EH4214	Human KYN(Kynurenine) ELISA Kit	7.813-500pmol/ml	ER1507	Rat E2(Estradiol) ELISA Kit	12.5-800pg/ml
EH4235	Human GAGs(Glycosaminoglycan) ELISA Kit	1.563-100ng/ml	ER0728	Rat ALP(Alkaline Phosphatase) ELISA Kit	0.625-40ng/ml
EM0010	Mouse MPO(Myeloperoxidase) ELISA Kit	0.156-10ng/ml	EM1503	Mouse PGE2(Prostaglandin E2) ELISA Kit	31.25-2000pg/ml
EM0109	Mouse IL-1β(Interleukin 1 Beta) ELISA Kit	12.5-800pg/ml	EM0843	Mouse ApoB100(Apolipoprotein B100) ELISA Kit	78.125-5000ng/ml
EM0121	Mouse IL-6(Interleukin 6) ELISA Kit	15.625-1000pg/ml	EM0514	Mouse Gal(Galanin peptides) ELISA Kit	31.25-2000pg/ml

## Hot Selling ELISA Kits

Product Code	Product Name	Range	Product Code	Product Name	Range
EH0012	Human Lipocalin-2 ELISA Kit	0.313-20ng/ml	EH4040	Human VDACC1 ELISA Kit	31.25-2000pg/ml
EH0017	Human MMP-2 ELISA Kit	0.5-32ng/ml	EH4047	Human DLL4(Delta-like protein 4) ELISA Kit	78.125-5000pg/ml
EH0018	Human MPO(Myeloperoxidase) ELISA Kit	1.563-100ng/ml	EH4136	Human cFN(Cellular Fibronectin) ELISA Kit	1.563-100ng/ml
EH0022	Human PLGF ELISA Kit	15.625-1000pg/ml	EH4140	Human calprotectin ELISA Kit	0.156-10ng/ml
EH0024	Human VEGFR1 ELISA Kit	125-8000pg/ml	EH4156	Human 8-iso-PGF2α(8-isoprostane) ELISA Kit	15.625-1000pg/ml
EH0025	Human A2M(Alpha-2-macroglobulin) ELISA Kit	78.125-5000ng/ml	EH4159	Human Visfatin C termina ELISA Kit	0.313-20ng/ml
EH0027	Human ACE2 ELISA Kit	0.391-25ng/ml	EH4161	Human anti-nuclear antibody ELISA Kit	1.563-100ng/ml
EH0028	Human Activin A ELISA Kit	62.5-4000pg/ml	EH4166	Human anti-NMDAR Antibody ELISA Kit	1.563-100ng/ml
EH0033	Human AR(Androgen receptor) ELISA Kit	0.313-20ng/ml	EH4191	Human ULK1(Unc-51 Like Kinase 1) ELISA Kit	15.625-1000pg/ml
EH0042	Human BAFF ELISA Kit	62.5-4000pg/ml	EH4243	Human SBDP 145 ELISA Kit	0.625-40ng/ml
EH0043	Human BDNF ELISA Kit	31.25-2000pg/ml	EH4245	Human SP(Substance P) ELISA Kit	7.813-500pg/ml
EH0051	Human CTSB(Cathepsin B) ELISA Kit	0.156-10ng/ml	EH4250	Human 3MH(3-methylhistidine) ELISA Kit	6.25-400nmol/ml
EH0092	Human CHEM(Chemerin) ELISA Kit	0.156-10ng/ml	EH4255	Human pro-BDNF ELISA Kit	15.625-1000pg/ml
EH0093	Human CHI3L1 ELISA Kit	62.5-4000pg/ml	EH4271	Human RF-IgM ELISA Kit	4.688-300IU/ml
EH0095	Human COMP ELISA Kit	0.625-40ng/ml	EH4274	Human CAMP ELISA Kit	0.469-30ng/ml
EH0102	Human CXCL10 ELISA Kit	7.813-500pg/ml	EH4278	Human SAKL ELISA Kit	93.75-6000pg/ml
EH0105	Human CXCL16 ELISA Kit	0.156-10ng/ml	EH4320	Human OVA(Ovalbumin) ELISA Kit	0.781-50ng/ml
EH0106	Human CXCL5 ELISA Kit	31.25-2000pg/ml	EH4328	Human IL-17 F(Interleukin 17F) ELISA Kit	15.625-1000pg/ml
EH0110	Human CST3(Cystatin C) ELISA Kit	0.313-20ng/ml	EH4343	Human SCUBE2 ELISA Kit	31.25-2000pg/ml
EH0125	Human ESM1 ELISA Kit	15.625-1000pg/ml	EH4349	Human Spexin ELISA Kit	7.813-500pg/ml
EH0128	Human FETUA(Fetuin A) ELISA Kit	0.781-50ng/ml	EH4702	Human Irisin ELISA Kit	1.563-100ng/ml
EH0130	Human FGF21 ELISA Kit	31.25-2000pg/ml	EH4820	Human CAF ELISA Kit	46.875-3000pg/ml
EH0134	Human FN1(Fibronectin) ELISA Kit	1.563-100ng/ml	EH4858	Human OxHDL ELISA Kit	1.563-100ng/ml
EH0140	Human FOLR1(Folate receptor alpha) ELISA Kit	15.625-1000pg/ml	EH7221	Human CD81(CD81 antigen) ELISA Kit	0.313-20ng/ml
EH0141	Human FRCTalkine ELISA Kit	125-8000pg/ml	EH7533	Human COLI(Collagen Type I) ELISA Kit	0.313-20ng/ml
EH0143	Human FUR(Furin) ELISA Kit	15.625-1000pg/ml	EM1522	Mouse PCT(Procalcitonin) ELISA Kit	7.813-500pg/ml
EH0145	Human GAL3(Galectin-3) ELISA Kit	0.156-10ng/ml	EM1532	Mouse esRAGE ELISA Kit	7.813-500pg/ml
EH0148	Human GAL9(Galectin-9) ELISA Kit	0.156-10ng/ml	EM1570	Mouse Mcpt1(Mast cell protease 1) ELISA Kit	78.125-5000pg/ml
EH0149	Human G-CSF ELISA Kit	39.063-2500pg/ml	EM1583	Mouse 4-HNE(4-Hydroxynonenal) ELISA Kit	31.25-2000pg/ml
EH0150	Human GDF15 ELISA Kit	23.438-1500pg/ml	EM1588	Mouse NLRP3 ELISA Kit	0.313-20ng/ml
EH0159	Human HGF ELISA Kit	125-8000pg/ml	EM1592	Mouse IL-8/cxcl15(Interleukin 8) ELISA Kit	15.625-1000pg/ml
EH0161	Human sICAM-1/CD54 ELISA Kit	0.781-50ng/ml	EM1607	Mouse ANA(Anti-nuclear Antibody) ELISA Kit	0.469-30ng/ml
EH0164	Human IFN-γ(Interferon gamma) ELISA Kit	15.625-1000pg/ml	EM1620	Mouse Calprotectin(Calprotectin) ELISA Kit	7.813-500ng/ml
EH0165	Human IGF-1 ELISA Kit	62.5-4000pg/ml	EM1710	Mouse Mcpt4(Mast cell protease 4) ELISA Kit	46.875-3000pg/ml
EH0166	Human IGF-2 ELISA Kit	62.5-4000pg/ml	EM1718	Mouse MASP2 ELISA Kit	0.781-50ng/ml
EH0169	Human IGFBP-3 ELISA Kit	0.391-25ng/ml	EM1850	Mouse Testosterone ELISA Kit	0.313-20ng/ml
EH0171	Human IGFBP-7 ELISA Kit	31.25-2000pg/ml	EM1853	Mouse FLG ELISA Kit	15.625-1000pg/ml
EH0173	Human IL-10(Interleukin 10) ELISA Kit	7.813-500pg/ml	EM1862	Mouse KYN(Kynurenine) ELISA Kit	7.813-500pmol/ml
EH0174	Human IL-11(Interleukin 11) ELISA Kit	31.25-2000pg/ml	ER1393	Rat TNF-α ELISA Kit	3.906-250pg/ml
EH0176	Human IL-12(Interleukin 12) ELISA Kit	15.625-1000pg/ml	ER0073	Rat Alb(Serum albumin) ELISA Kit	7.813-500ng/ml
EH0177	Human IL-15(Interleukin 15) ELISA Kit	15.625-1000pg/ml	ER1487	Rat 8-OHdG ELISA Kit	1.563-100ng/ml
EH0184	Human IL-1α(Interleukin 1 Alpha) ELISA Kit	3.906-250pg/ml	ER1618	Rat CB1R/CNR1 ELISA Kit	78.125-5000pg/ml
EH0185	Human IL-1β(Interleukin 1 Beta) ELISA Kit	3.906-250pg/ml	ER1630	Rat IL-22(Interleukin 22) ELISA Kit	15.625-1000pg/ml
EH0198	Human IL-33(Interleukin-33) ELISA Kit	15.625-1000pg/ml	ER1636	Rat Calprotectin(Calprotectin) ELISA Kit	7.813-500ng/ml
EH0201	Human IL-6 ELISInterleukin-6 receptor sA Kit	4.688-300pg/ml	ER1637	Rat ANG II(Angiotensin II) ELISA Kit	31.25-2000pg/ml
EH0202	Human IL-6Rα(ubunit alpha) ELISA Kit	31.25-2000pg/ml	ER1707	Rat GABA ELISA Kit	31.25-2000pg/ml
EH0203	Human IL-7(Interleukin 7) ELISA Kit	7.813-500pg/ml	ER1800	Rat PGE2(Prostaglandin E2) ELISA Kit	31.25-2000pg/ml
EH0205	Human IL-8(Interleukin 8) ELISA Kit	31.25-2000pg/ml	EU0192	PAF(Platelet Activating Factor) ELISA Kit	0.156-10ng/ml
EH0209	Human SERPINAA4(Kallistatin) ELISA Kit	15.625-1000ng/ml	EU0197	HCV(Homocysteine) ELISA Kit	7.813-500pmol/ml
EH0216	Human LEP(Leptin) ELISA Kit	31.25-2000pg/ml	EU0198	AP(Allopregnanolone) ELISA Kit	1.563-100ng/ml
EH0222	Human MCP-1 ELISA Kit	15.625-1000pg/ml	EU0199	MT(Melatonin) ELISA Kit	7.813-500pg/ml
EH0232	Human MMP-1 ELISA Kit	78.125-5000pg/ml	EU0382	HIS(Histamine) ELISA Kit	1.563-100ng/ml
EH0237	Human MMP-8 ELISA Kit	0.156-10ng/ml	EU0187	4-HNE(4-Hydroxynonenal) ELISA Kit	31.25-2000pg/ml
EH0238	Human MMP-9 ELISA Kit	0.313-20ng/ml	EU0188	KYN(Kynurenine) ELISA Kit	7.813-500pmol/ml
EH0255	Human POSTN/OSF2(Periostin) ELISA Kit	0.156-10ng/ml	EU3118	PD1(Protectin D1) ELISA Kit	0.156-10ng/ml
EH0265	Human Rantes ELISA Kit	31.25-2000pg/ml	EH4869	Human Cd19 ELISA Kit	0.156-10ng/ml
EH0266	Human RBP4 ELISA Kit	1.563-100ng/ml	ER1720	Rat T3(Triiodothyronine) ELISA Kit	0.156-10ng/ml
EH0269	Human RETN(Resistin) ELISA Kit	31.25-2000pg/ml	ER1721	Rat T4(Thyroxine) ELISA Kit	0.703-45ng/ml
EH0278	Human SDC1(Syndecan-1) ELISA Kit	0.156-10ng/ml	EH12498	Human SOS1 ELISA Kit	0.156-10ng/ml
EH0288	Human TGF-β2 ELISA Kit	15.625-1000pg/ml	EH15230	Human RIPK3 ELISA Kit	0.313-20ng/ml

## Hot Selling Antibodies

Product Code	Product Name	Reactivity	Application
FNab00129	anti- ACVR1B antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00839	anti- BCL2 antibody	Human, Mouse, Rat	ELISA, WB, IHC, IF
FNab00859	anti- BDNF antibody	Human, Rat, Mouse	ELISA, WB, IHC, IF
FNab00869	anti- beta actin antibody	Human, Mouse, Rat	ELISA, IF, WB, IHC
FNab00917	anti- BMP2 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab01469	anti- CD34 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab02088	anti- CX3CL1 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab02158	anti- CYP2D6 antibody	Human	ELISA, WB
FNab02941	anti- FADD antibody	Human, Mouse, Rat	ELISA, IHC, WB, IP
FNab03342	anti- GAPDH antibody	Human, Mouse, Rat	ELISA, IHC, WB
FNab03343	anti- GAPDH antibody	Human, Mouse, Rat, Zebrafish	ELISA, IHC, WB, IF, FC, IP
FNab03924	anti- HMGB1 antibody	Human, Mouse, Rat	ELISA, WB, IHC, IF
FNab04279	anti- IL-4 antibody	Human, Rat	ELISA, WB, IHC
FNab04325	anti- iNOS antibody	Human, Mouse, Rat	ELISA, WB, IHC, IF
FNab04551	anti- KIF15 antibody	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP
FNab05026	anti- MAT1A antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab05100	anti- MEF2A antibody	Human, Mouse, Rat	ELISA, WB
FNab05238	anti- MMP2 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab05247	anti- MMP9 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab05579	anti- NCAPH antibody	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP
FNab05759	anti- NLRX1 antibody	Human, Mouse, Rat	ELISA, WB, IF, IHC, IP
FNab05764	anti- NME3 antibody	Human, Mouse, Rat	ELISA, WB
FNab05850	anti- NRBF2 antibody	Human, Rat	ELISA, WB, IHC, IF
FNab06089	anti- p65/RELA antibody	Human, Mouse, Rat, Pig	ELISA, IHC, IF, IP, WB, FC
FNab06712	anti- PPP2CA antibody	Human, Mouse, Rat	ELISA, IHC, WB, IF
FNab06912	anti- PTEN antibody	Human, Mouse	ELISA, IHC, WB, IF
FNab07515	anti- RtcB-Specific antibody	Human, Mouse	ELISA, WB, IF, IP
FNab08124	anti- SOX17 antibody	Human, Mouse, Rat	ELISA, WB
FNab08216	anti- SRC antibody	Human, Mouse, Rat, Hamster	ELISA, IHC, IP, FC, IF, WB
FNab08310	anti- Stathmin antibody	Human, Mouse	ELISA, WB, IHC, IF
FNab08334	anti- STK33 antibody	Human, Mouse, Rat	ELISA, WB, IHC, IP
FNab08821	anti- TNF-alpha antibody	Human, Mouse, Rat	ELISA, IHC, WB
FNab09285	anti- URAT1 antibody	Human, Mouse, Rat	ELISA, WB
FNab09380	anti- VCAM-1 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab09559	anti- YAP1 antibody	Human, Mouse	ELISA, WB, IHC, IF, IP
FNab09774	anti- PR antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab09797	anti- Asprosin antibody	Human, Mouse	ELISA, IHC, WB
FNab09805	anti- P Glycoprotein antibody	Human, Mouse	ELISA, WB, IP, IHC
FNab09875	anti- p65/RELA antibody	Human, Mouse, Rat	ELISA, IHC, IF, WB
FNab10013	anti- Caspase 3 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab10034	anti- CX3CR1 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab10218	anti- HMGB1 antibody	Human	ELISA, WB
FNab10269	anti- OAT1 antibody	Human, Mouse	ELISA, WB
FNab10270	anti- GLUT9 antibody	Human, Mouse	ELISA, WB, IHC
FNab00129	anti- ACVR1B antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab00839	anti- BCL2 antibody	Human, Mouse, Rat	ELISA, WB, IHC, IF
FNab00859	anti- BDNF antibody	Human, Rat, Mouse	ELISA, WB, IHC, IF
FNab00869	anti- beta actin antibody	Human, Mouse, Rat	ELISA, IF, WB, IHC
FNab00917	anti- BMP2 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab01469	anti- CD34 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab02088	anti- CX3CL1 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab02158	anti- CYP2D6 antibody	Human	ELISA, WB
FNab02941	anti- FADD antibody	Human, Mouse, Rat	ELISA, IHC, WB, IP
FNab03342	anti- GAPDH antibody	Human, Mouse, Rat	ELISA, IHC, WB
FNab03343	anti- GAPDH antibody	Human, Mouse, Rat, Zebrafish	ELISA, IHC, WB, IF, FC, IP
FNab03924	anti- HMGB1 antibody	Human, Mouse, Rat	ELISA, WB, IHC, IF
FNab04279	anti- IL-4 antibody	Human, Rat	ELISA, WB, IHC
FNab04325	anti- iNOS antibody	Human, Mouse, Rat	ELISA, WB, IHC, IF
FNab04551	anti- KIF15 antibody	Human, Mouse, Rat	ELISA, WB, IHC, IF, IP
FNab05026	anti- MAT1A antibody	Human, Mouse, Rat	ELISA, WB, IHC

## Hot Selling Antibodies

Product Code	Product Name	Reactivity	Application
FNab09543	anti- XK antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab09559	anti- YAP1 antibody	Human, Mouse	ELISA, WB, IHC, IF, IP
FNab09572	anti- YTHDF1 antibody	Human, Mouse, Rat	ELISA, WB, IHC, IP
FNab09751	anti- ZO-1 antibody	Human, Mouse, Rat	ELISA, WB, IHC, FC, IF, IP
FNab09780	anti- CD45 antibody	Human, Pig	ELISA, WB, IHC, FC
FNab09788	anti- KI67 antibody	Human, Rat, Mouse	ELISA, WB, IHC, IF
FNab09789	anti- Elastin antibody	Human, Rat, Mouse	ELISA, IHC, IP, WB
FNab09797	anti- Asprosin antibody	Human, Mouse	ELISA, IHC, WB
FNab09800	anti- Vinculin antibody	Human, Rat, Mouse	ELISA, FC, IF, IHC, WB
FNab09812	anti- CD206 antibody	Human, Mouse, Rat	ELISA, IP, WB
FNab09837	anti- TLR4 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab09851	anti- TGFβ1 antibody	Human, Mouse, Rat	ELISA, WB
FNab09856	Anti-Tyrosinase antibody	Human, Mouse, Rat	ELISA, WB, IHC, IF
FNab09864	anti- CXCL8/IL8 antibody	Human, Mouse	ELISA, IHC, WB
FNab09870	anti- TH antibody	Human, Mouse, Rat	ELISA, WB, IHC, IF
FNab09872	anti- Beta Tubulin antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab09875	anti- p65/RELA antibody	Human, Mouse, Rat	ELISA, IHC, IF, WB
FNab09891	anti- alpha-synuclein antibody	Human, Mouse, Rat	ELISA, IHC, WB
FNab09509	anti- WFS1 antibody	Human, Mouse	ELISA, IHC, IF, WB
FNab09524	anti- WRN antibody	Human	ELISA, WB
FNab09936	anti- IGF1 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab09938	anti- CASP3 antibody	Human, Mouse, Rat	ELISA, IHC, WB, IF
FNab09967	anti- CLDN5 antibody	Human, Mouse, Rat	ELISA, WB, IHC, IF
FNab09972	anti- PDGFRβ antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab09983	anti- CBPN antibody	Human	ELISA, WB, IHC
FNab10013	anti- Caspase 3 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab10014	anti- BDNF antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab10015	anti- Caspase 1 antibody	Human, Mouse	ELISA, WB, IHC
FNab10034	anti- CX3CR1 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab10038	anti- KLF4 antibody	Human, Rat, Mouse	ELISA, WB, IHC
FNab10088	anti- Parkin antibody	Human, Rat	ELISA, WB, IHC
FNab10113	anti- MTNR1B antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab10120	anti- NLRP3 antibody	Human, Mouse, Rat	ELISA, WB, IHC, IF
FNab10121	anti- CD200 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab10218	anti- HMGB1 antibody	Human	ELISA, WB
FNab10227	anti- akt1 antibody	Human, rat	ELISA, WB, IHC
FNab10245	anti- FOXO3A antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab10250	anti- HSP60 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab10258	anti- PTEN antibody	Human, Rat	ELISA, WB
FNab10331	anti- p38 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab10334	anti- NFKB1 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab10340	anti- HIF1a antibody	Human, Mouse, Rat	ELISA, IHC, WB
FNab10342	anti- c-MYC antibody	Human, Mouse, Rat	ELISA, WB
FNab10407	anti- COX2 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab10408	anti- CD34 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab10410	anti- BAX antibody	Human	ELISA, WB
FNab10412	anti- BACE1 antibody	Human, Rat	ELISA, WB, IHC
FNab10429	anti- CD68 antibody	Human, Mouse, Rat	ELISA, IHC, WB
FNab10437	anti- CXCR4 antibody	Human, Rat	ELISA, WB, IHC
FNab10441	anti- γ-H2AX antibody	Human, Mouse, Rat	WB, IHC, IF
FNab10448	anti- CD81 antibody	Human, Mouse, Rat	ELISA, IHC, WB
FNab10452	anti- GPX4 antibody	Human, Rat	ELISA, WB, IHC
FNab10467	anti- Cyclin D1 antibody	Human, Mouse	ELISA, WB, IHC
FNab10472	Anti-NCR1 antibody	Human, Mouse	ELISA, WB, IHC, IF
FNab10476	Anti-CD137 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab10536	anti- GRIN1 antibody	Human, Mouse, Rat	ELISA, WB, IHC
FNab10699	Goat anti- Monkey IgG(gamma)	monkey	ELISA, WB
FNab10710	anti- HMGB1 antibody	Human, Mouse, Rat	ELISA, WB, IHC, IF
FNab10724	anti- PIK3R1 antibody	Human, Mouse	ELISA, WB
FNab10728	anti- TK1 antibody	Human	ELISA, WB, IF, IHC

## Hot Selling Secondary Antibodies

Product Code	Product Name	Application
Alexa Fluor 405 conjugated Affinipure Goat Anti-Mouse IgG(H+L)	FNSA-0094	IF, FC, IHC
Alexa Fluor 488 conjugated Affinipure Goat Anti-Mouse IgG(H+L)	FNSA-0055	IF, FC
Alexa Fluor 488 conjugated Affinipure Goat Anti-Rabbit IgG(H+L)	FNSA-0056	IF, FC
Alexa Fluor 555 conjugated Affinipure Goat Anti-Rabbit IgG(H+L)	FNSA-0095	IF, FC
Alexa Fluor 594 conjugated Affinipure Goat Anti-Mouse IgG(H+L)	FNSA-0060	IF, FC
Alexa Fluor 594 conjugated Affinipure Goat Anti-Rabbit IgG(H+L)	FNSA-0061	IF, FC
AMCA conjugated Affinipure Goat Anti-Mouse IgG (H+L)	FNSA-0082	IF, FC
AMCA conjugated Affinipure Goat Anti-Rabbit IgG (H+L)	FNSA-0083	IF, FC
AP conjugated Affinipure Goat Anti- Rat IgG (H+L)	FNSA-0028	ELISA, WB
AP conjugated Affinipure Goat Anti-Mouse IgG (H+L)	FNSA-0020	ELISA, WB
Biotin conjugated 2019 nCOV N protein antibody(1G5)	FNSA-0115	ELISA, WB, IHC
Biotin conjugated 2019 nCOV N protein antibody(3C6)	FNSA-0116	ELISA, WB, IHC
Biotin conjugated Affinipure Goat Anti-Human IgG (H+L)	FNSA-0044	ELISA, WB, IHC/ICC
Biotin conjugated Affinipure Goat Anti-Mouse IgG(H+L)	FNSA-0040	ELISA, WB, IHC/ICC
Biotin conjugated Affinipure Goat Anti-Rabbit IgG (H+L)	FNSA-0041	ELISA, WB, IHC/ICC
Cy3 conjugated Affinipure Goat Anti-Mouse IgG(H+L)	FNSA-0077	IF, Flow cyt
Cy3 conjugated Affinipure Goat Anti-Rabbit IgG (H+L)	FNSA-0078	IF, Flow cyt
FITC conjugated Affinipure Goat Anti Human IgA (H+L)	FNSA-0096	IF, Flow cyt
FITC conjugated Affinipure Goat Anti Human IgG (H+L)	FNSA-0038	IF, Flow cyt
FITC conjugated Affinipure Goat Anti Human IgM (H+L)	FNSA-0097	IF, Flow cyt
FITC conjugated Affinipure Goat Anti Mouse IgG (H+L)	FNSA-0030	IF, Flow cyt
FITC conjugated Affinipure Goat Anti Rabbit IgG (H+L)	FNSA-0031	IF, Flow cyt
HRP conjugated Goat anti Mouse IgG1	FNSA-0086	ELISA, WB
HRP conjugated Goat anti Mouse IgG2a	FNSA-0087	ELISA, WB
HRP conjugated Goat anti Mouse IgG2b	FNSA-0088	ELISA, WB
HRP conjugated Goat anti Mouse IgG2c	FNSA-0089	ELISA, WB
HRP conjugated Goat anti Mouse IgG3	FNSA-0090	ELISA, WB
HRP conjugated Goat anti Mouse IgM	FNSA-0091	ELISA, WB
poly-HRP conjugated Chicken anti rabbit IgG (H+L)	FNSA-0118	ELISA, WB
poly-HRP conjugated Mouse anti- human IgG	FNSA-0114	ELISA
poly-HRP conjugated Mouse anti- human IgM	FNSA-0113	ELISA
Rhodamine (TRITC) conjugated Affinipure Goat Anti-Human IgG (H+L)	FNSA-0070	IF, Flow cyt
Rhodamine (TRITC) conjugated Affinipure Goat Anti-Mouse IgG(H+L)	FNSA-0065	IF, Flow cyt
Rhodamine (TRITC) conjugated Affinipure Goat Anti-Rabbit IgG (H+L)	FNSA-0066	IF, Flow cyt

## Hot Selling Proteins

Product Code	Product Name	Species	Source	Tags	Product Code	Product Name	Species	Source	Tags
P0010	Recombinant Human PCT	Human	E.Coli	His tag	P0135	Recombinant Human Casp3	Human	E.Coli	His tag
P0042	Recombinant Human Insulin	Human	E.Coli	His tag	P0136	Recombinant Bovine Uter	Bovine	E.Coli	His tag
P0055	Recombinant Human Asprosin	Human	E.Coli	His tag	P0137	Recombinant Human Uter	Human	E.Coli	His tag
P0067	Recombinant Human $\beta$ -Actin	Human	E.Coli	His tag	P0138	Recombinant Hamster Uter	Hamster	E.Coli	His tag
P0108	Recombinant Human HSPA5	Human	E.Coli	His tag	P0139	Recombinant Mouse Uter	Mouse	E.Coli	His tag
P0109	Recombinant Human CXCL1	Human	E.Coli	His tag	P0140	Recombinant Rat Uter	Rat	E.Coli	His tag
P0113	Recombinant EGFP protein	Jellyfish	E.Coli	His tag	P0141	Recombinant Bovine PAI-1	Bovine	E.Coli	His tag
P0116	Recombinant Human GAPDH	Human	E.Coli	His tag	P0142	Recombinant Human PAI-1	Human	E.Coli	His tag
P0117	Recombinant Human YY1	Human	E.Coli	His tag	P0143	Recombinant Mouse PAI 1	Mouse	E.Coli	His tag
P0119	Recombinant Human trx-1	Human	E.Coli	His tag	P0144	Recombinant Pig PAI-1	Pig	E.Coli	His tag
P0120	Recombinant Mouse txn	Mouse	E.Coli	His tag	P0145	Recombinant Rat PAI-1	Rat	E.Coli	His tag
P0122	Recombinant Human FABP2	Human	E.Coli	His tag	P0146	Recombinant Bovine IFNG	Bovine	E.Coli	His tag
P0124	Recombinant Human FABP4	Human	E.Coli	His tag	P0147	Recombinant Chicken IFNG	Chicken	E.Coli	His tag
P0125	Recombinant Human FABP5	Human	E.Coli	His tag	P0148	Recombinant Dog IFNG	Dog	E.Coli	His tag
P0126	Recombinant Human FABP6	Human	E.Coli	His tag	P0149	Recombinant Human IFNG	Human	E.Coli	His tag
P0127	Recombinant Human FABP7	Human	E.Coli	His tag	P0150	Recombinant Mouse IFNG	Mouse	E.Coli	His tag
P0128	Recombinant Human FABP8	Human	E.Coli	His tag	P0151	Recombinant Pig IFNG	Pig	E.Coli	His tag
P0129	Recombinant Human FABP9	Human	E.Coli	His tag	P0152	Recombinant Rat IFNG	Rat	E.Coli	His tag
P0130	Recombinant Mouse FABP4	Mouse	E.Coli	His tag	P0153	Recombinant Bovine TNFA	Bovine	E.Coli	His tag
P0132	Recombinant Human CYGB	Human	E.Coli	His tag	P0154	Recombinant Dog TNFA	Dog	E.Coli	His tag

## Selected Publications

Currently, nearly 1,000 academic publications citing FineTest products are published in well-known journals. Now, hundreds of newly cited papers are increasing every year.



Product Code	Product Name	IF	Journal Title	Article Title
EH4233	Human PGE2(Prostaglandin E2) ELISA Kit	31.398	Cell	Acetylation Blocks cGAS Activity and Inhibits Self-DNA-Induced Autoimmunity
EH2484	Human RHOA(Transforming protein RhoA)	16.836	Advanced Functional Materials	Tumor Microenvironment - Activatable Nanoenzymes for Mechanical Remodeling of Extracellular Matrix and Enhanced Tumor Chemotherapy
EH2259	Human PRKCSH(Glucosidase 2 subunit beta) ELISA Kit	16.836	Nanotoday	Nano-scarifiers for blood biomarker discovery in ovarian carcinoma
EH9300	Human IFITM2 ELISA Kit	14.079	hepatology	Exosomal IFITM2 transmitted to DCs inhibits IFN $\alpha$ pathway activation and blocks anti - HBV efficacy of exogenous IFN $\alpha$
FNab05465	anti-Myelin basic protein-antibody	13.334	Acs Nano	Flexible and Highly Biocompatible Nanober- Based Electrodes for NeuralS urface Interfacing
EH3057	Human FG(Fibrinogen) ELISA Kit	11.337	Brain	Systemic infection exacerbates cerebrovascular dysfunction in Alzheimer's disease
EH1580	Human ACLY(ATP-citRate synthase) ELISA Kit	9.727	Cancer Research	Long non-coding RNA TINCR-mediated regulation of acetyl-CoA metabolism promotes nasopharyngeal carcinoma progression and chemoresistance
EM0382	Mouse Hmgb1(High mobility group protein B1) ELISA Kit	8.728	Journal for ImmunoTherapy of Cancer	Cell death induced by cytotoxic CD8+ T cells is immunogenic and primes caspase-3-dependent spread immunity against endogenous tumor antigens.
EH1845	Human PTMA(Prothymosin alpha) ELISA Kit	8.578	The Journal of Heart and Lung Transplantation	prothymosin $\alpha$ gene transduction attenuated cardiopulmonary remodeling and mortality in flow-induced pulmonary hypertension rat model
EH8929	Human GRASP65/GORASP1 ELISA Kit	8.456	ACS Applied Materials & Interfaces	supramolecular insights into domino effects of ag@zno-induced oxidative stress in melanoma cancer cells
EH4235	Human GAGs(Glycosaminoglycan) ELISA Kit	8.308	CHEST	a cluster analysis of bronchiectasis patients based on the airway immune profile
EH0201	Human IL-6(Interleukin 6) ELISA Kit	8.084	World Journal of Pharmaceutical Research	comparative evaluation of povidone iodine and aloe vera on gingival health: a randomized controlled trial
EH1534	Human GPLD1(Phosphatidylinositol-glycan-specific phospholipase D) ELISA Kit	7.919	Clinical and Translational Medicine	Reiterative infusions of MSCs improve pediatric osteogenesis imperfecta eliciting a pro-osteogenic paracrine response: TERCELOI clinical trial
EM1324	Mouse PYGL (Glycogen Phosphorylase, Liver) ELISA Kit	7.815	Cell Reports	Glycogen Dynamics Drives Lipid Droplet Biogenesis during Brown Adipocyte Differentiation
EM0260	Mouse INS(Insulin) ELISA Kit	7.65	Journal of Hazardous Materials	Arsenic induces pancreatic dysfunction and ferroptosis via mitochondrial ROS-autophagy-lysosomal pathway
EM0088	Mouse GDNF(Glial Cell Line Derived Neurotrophic Factor) ELISA Kit	7.242	Acta Biomaterialia	Biomaterial based strategies to reconstruct the nigrostriatal pathway in organotypic slice co-cultures
EM1400	Mouse TFR(Transferrin Receptor) ELISA Kit	7.116	Haematologica	Cell-specific expression of Hfe determines the outcome of Salmonella enterica serovar Typhimurium infection in mice