

FC Auxiliary Reagents

Product Category	Cat.No	Product Name	Size
Stimulating and Blocking	K083	Cell Stimulation and Protein Transport Inhibitor Kit	50T/100T/200T
Fixation/Permeabilization	K081	Foxp3/Transcription Factor Staining Kit	20T
Fixation/Permeabilization	K082R	Intracellular Fixation/Permeabilization Buffer Kit (Ready To Use)	50T/100T/500T
FcR Blocking	FNab30106	Purified Anti-Mouse CD16/32 Antibody	100µg
FcR Blocking	FNab10997	Purified Anti-Human CD16 Antibody	100µg
Dead Double Staining	FNCK022	Calcein AM/PI Live/Dead Double Staining Kit	100T/500T
Dead Double Staining	K121	FineTest® Red 780(Viability Dye)	100T

Isotype Control Antibodies

Cat.No	Product Name	Host	Marker	Size
FITC-30073	FITC Rat IgG2a,κ Isotype Control(2A3)	Rat	FITC	20T/100T
PE-30073	PE Rat IgG2a,κ Isotype Control(2A3)	Rat	PE	20T/100T
APC-30073	APC Rat IgG2a,κ Isotype Control(2A3)	Rat	APC	20T/100T
FITC-30075	FITC Rat IgG1,κ Isotype Control(HRPN)	Rat	FITC	20T/100T
PE-30075	PE Rat IgG1,κ Isotype Control(HRPN)	Rat	PE	20T/100T
APC-30075	APC Rat IgG1,κ Isotype Control(HRPN)	Rat	APC	20T/100T
FITC-30076	FITC Mouse IgG1,κ Isotype Control(MOPC-21)	Mouse	FITC	20T/100T
PE-30076	PE Mouse IgG1,κ Isotype Control(MOPC-21)	Mouse	PE	20T/100T
APC-30076	APC Mouse IgG1,κ Isotype Control(MOPC-21)	Mouse	APC	20T/100T
PCP-30076	PerCP Mouse IgG1,κ Isotype Control(MOPC-21)	Mouse	PerCP	20T/100T
FITC-30077	FITC Mouse IgG2a,κ Isotype Control(C1.18.4)	Mouse	FITC	20T/100T
PE-30077	PE Mouse IgG2a,κ Isotype Control(C1.18.4)	Mouse	PE	20T/100T
APC-30077	APC Mouse IgG2a,κ Isotype Control(C1.18.4)	Mouse	APC	20T/100T
FITC-30078	FITC Mouse IgG2b,κ Isotype Control(MPC-11)	Mouse	FITC	20T/100T
PE-30078	PE Mouse IgG2b,κ Isotype Control(MPC-11)	Mouse	PE	20T/100T
APC-30078	APC Mouse IgG2b,κ Isotype Control(MPC-11)	Mouse	APC	20T/100T



Wuhan Fine Biotech Co., Ltd

1-2 Floor, BLD 22, Optics Valley Biopharmaceutical Accelerator, No.388  
Gaoxin 2nd Road, East Lake High Tech Development Zone, Wuhan,  
430074, Hubei, China

[www.fn-test.com](http://www.fn-test.com)

[fine@fn-test.com](mailto:fine@fn-test.com)



Wuhan Fine Biotech Co., Ltd.



Flow Cytometry Antibodies





● Applications: Widely used in the fields of immunology, medicine and biology etc, such as cell phenotype identification, cell sorting, cytokine measurement, cell proliferation, cell apoptosis, cell cycle, cell killing etc.

★ Flow Cytometry Antibodies

An efficient assistant for immune cell identification and a convenient tool for cell status identification.

FineTest® provides better experience with high-quality flow cytometry antibodies. The products cover human, mouse, rat and so on, with rich targets. FineTest® offers you an integrative flow cytometry solution including free colour matching, assay guidance and data analysis services etc.



Advantages

01 High Quality

Derived from classically cloned and excellently labelling technology

02 Various Dyes


high-performance in multicolour flow cytometry

03 Complete Solution



free colour matching, assay guidance and data analysis services

Multiple Fluorescein Labeled

More than 10 kinds of fluoresceins, which meet the requirements of color matching in multi-color flow experiments.



**FITC-30001**  
**FITC Anti-Mouse CD4 (GK1.5)**  
**Results :**



Application Of Cell Grouping

T/B/NK  
Cells

DC  
Cells

Th1  
Cells

Th2  
Cells

Th17  
Cells

Treg  
Cells

Macrophage  
cells

Stem  
cells

Flow Cytometry Colour Matching Recommendation

Species	Cell types	Recommendation	Cat.No
Human	T/B/NK Cells	CD45-PerCP	PCP-30039
		CD3-FITC	FITC-30004
		CD16-PE	PE-30061
		CD56-PE	PE-30008
		CD19-APC	APC-30066
Human	Th1/Th2 Cells	CD3-PerCP/Cyanine5.5	PCP55-30004
		CD4-FITC	FITC-30005
		IFN-γ-PE	PE-30053
		IL4-APC	APC-30043
Mouse	Th1/Th2 Cells	CD3-PerCP/Cyanine5.5	PCP55-30002
		CD4-FITC	FITC-30128
		IFN-γ-PE	PE-30074
		IL4-APC	APC-30026
Human	Treg Cells	CD4-FITC	FITC-30005
		CD25-PE	PE-30035
		CD3-PerCP-Cy5.5	PCP55-30004
		CD127-FineTest®647	F647-30033
Mouse	Treg Cells	CD4-FITC	FITC-30128
		CD25-APC	APC-30017
		FOXP3-PE	PE-30111

FC Auxiliary Reagents

Product Category	Cat.No	Product Name	Size
Erythrocyte Lysis Buffer	K080	10x ACK Lysis Buffer(without fixative)	100mL/200mL/500mL
Buffer	K079	Cell Staining Buffer	100mL/200mL/500mL
Nuclear Staining	K116	PI Reagent	50T/100T/200T/500T
Nuclear Staining	K115	7-AAD Reagent	50T/100T/200T/500T
Nuclear Staining	K117	DAPI Reagent	50T/100T/200T/500T