

➤ Anti-PEG Antibody ELISA Kits

FineTest® has developed a series of ELISA Kits for studying of the immune response induced by PEG-modified proteins or drugs. The anti-PEG IgM and IgG antibodies produced in human, monkey, rat and rabbit were quantitatively tested.

Cat.No	Product Name	Target	Sensitivity	Range
EH5111	Human anti-PEG IgG ELISA Kit	PEG	0.375ng/ml	0.625-40ng/ml
EH5112	Human anti-PEG IgM ELISA Kit	PEG	0.188ng/ml	0.313-20ng/ml
EM2149	Mouse anti-PEG IgG ELISA Kit	PEG	0.469ng/ml	0.781-50ng/ml
EM2150	Mouse anti-PEG IgM ELISA Kit	PEG	0.188ng/ml	0.313-20ng/ml
ER2094	Rat anti-PEG IgG ELISA Kit	PEG	0.75ng/ml	1.25-80ng/ml
ER2095	Rat anti-PEG IgM ELISA Kit	PEG	0.094ng/ml	0.156-10ng/ml
ERB0231	Rabbit anti-PEG IgG ELISA Kit	PEG	0.375ng/ml	0.625-40ng/ml
ERB0232	Rabbit anti-PEG IgM ELISA Kit	PEG	0.094ng/ml	0.156-10ng/ml
EMK0221	Monkey anti-PEG IgG ELISA Kit	PEG	0.375ng/ml	0.625-40ng/ml
EMK0222	Monkey anti-PEG IgM ELISA Kit	PEG	0.188ng/ml	0.313-20ng/ml

➤ Host Cell Proteins(HCPs)ELISA Kits

HCPs residue must be strictly controlled during the production of antibody drugs. FineTest® has developed the following HCPs ELISA Kits to support your research.

Cat.No	Product Name	Sensitivity	Range
EH4435-CM	HEK 293 HCPs(HEK 293 host cell proteins) ELISA Kit	46.875ng/ml	78.125-5000ng/ml
EMK0179	VERO HCPs(VERO host cell proteins) ELISA Kit	0.188ng/ml	0.313-20ng/ml
EU2624	E.coli HCPs(E.coli host cell proteins) ELISA Kit	1.875ng/ml	3.125-200ng/ml
EU2614	CHO HCPs ELISA Kit	0.188ng/ml	0.313-20ng/ml
EU3135	OVA(Ovalbumin) ELISA Kit	0.469ng/ml	0.781-50ng/ml
EU3133	PEA(Pseudomonas Exotoxin A) ELISA Kit	18.75pg/ml	31.25-2000pg/ml
EU3126	LPS(Lipopolysaccharides) ELISA Kit	0.188ug/ml	0.313-20ug/ml
EP0253	Pig TRY (Trypsin) ELISA Kit	0.938ng/ml	1.563-100ng/ml
EB0148	Bovine serum albumin(BSA) ELISA Kit	0.094ug/ml	0.156-10ug/ml
EU2662	Protein A ELISA Kit	46.875pg/ml	78.125-5000pg/ml
EU2607	FLAG-Tag (DYKDDDDK-Tag Protein) ELISA Kit	9.375ng/ml	15.625-1000ng/ml
EH4983	Human anti-cell free DNA ELISA Kit	0.938ng/ml	1.563-100ng/ml
EH4955	Human DS-DNA-IgG ELISA Kit	Qualitative	Qualitative



Wuhan Fine Biotech Co., Ltd

1-2 Floor, BLD 22, Optics Valley Biopharmaceutical Accelerator, No.388
Gaoxin 2nd Road, East Lake High Tech Development Zone, Wuhan,
430074, Hubei, China

www.fn-test.com

fine@fn-test.com



ELISA Kits Antibody Drugs



➤ ELISA Kits for therapeutic mAbs and ADAs

FineTest® ELISA kits are used for the detection and quantification of therapeutic monoclonal antibodies such as Bevacizumab® and biosimilars including rituximab® and adalimumab®, as well as the immunogenicity detection of Anti-Drug Antibodies (ADAs).

Advantages

High specificity

No cross-reactivity with non-target human serum protein or immunoglobulin

High sensitivity

High sensitivity and wide detection range

High accuracy

Low intra and inter assay variations

Monoclonal antibody drugs are exogenous proteins. Most of them have been modified via humanized heavychain, because Fab region is still immunogenic. Allergic and immune reactions often occur in large quantities or long-term injections. FineTest® has developed some ELISA Kits to help scientific researchers and medical workers to determine the metabolism of drugs during drug treatment.

Drug	Cat.No	Product Name	Sensitivity	Range
Adalimumab	EH5106	Human Adalimumab ELISA Kit	0.938ng/ml	1.563-100ng/ml
	EH5107	Human Anti-Adalimumab antibody ELISA kit	0.938ng/ml	1.563-100ng/ml
	EH5126	Human Anti-Adalimumab antibody ELISA kit	Qualitative	Qualitative
Aflibercept	EH5134	Human Aflibercept ELISA Kit	0.469ng/ml	0.781-50ng/ml
	EH5135	Human Anti-Aflibercept ELISA Kit	0.188ng/ml	0.313-20ng/ml
Atezolizumab	EU2709	Anti-Atezolizumab antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
	EU2713	Anti-Atezolizumab Neutralizing antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
Belimumab	EU7006	Belimumab ELISA kit	0.469ng/ml	0.781-50ng/ml
Benralizumab	EU2707	Anti-Benralizumab antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
	EU2708	Anti-Benralizumab Neutralizing antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
Bevacizumab	EU2594	Bevacizumab (Avastin) ELISA Kit	0.188ng/ml	0.313-20ng/ml
	EH5110	Human Anti-Bevacizumab(Avastin) antibody ELISA Kit	0.188ng/ml	0.313-20ng/ml
Burosumab	EU2704	Burosumab ELISA Kit	0.469ng/ml	0.781-50ng/ml
	EU2705	Anti-Burosumab antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
Daratumumab	EU6013	Anti-Daratumumab antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
Datopotamab	EU6985	Datopotamab deruxtecan (anti-TROP2) ELISA Kit	0.469ng/ml	0.781-50ng/ml
Denosumab	EH4515	Human Denosumab ELISA Kit	0.069ng/ml	0.137-100ng/ml
	EH5125	Human Anti-Denosumab antibody ELISA kit	0.938ng/ml	1.563-100ng/ml
	EH5127	Human Anti-Denosumab antibody ELISA kit	Qualitative	Qualitative
Donanemab	EH5289	Human Anti-Donanemab antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
Dupilumab	EH5290	Human Dupilumab ELISA Kit	0.938ng/ml	1.563-100ng/ml
	EH5291	Human Anti-Dupilumab antibody ELISA Kit	0.469ng/ml	0.781-50ng/ml
Elotuzumab	EU6729	Elotuzumab ELISA kit	0.469ng/ml	0.781-50ng/ml
Enfortumab	EU2636	Enfortumab (Anti-NECTIN4) ELISA Kit	0.188ng/ml	0.313-20ng/ml
	EH5294	Human Anti-Enfortumab antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
Faricimab	EU2639	Faricimab (Human VEGF/ang-2 antibody) ELISA Kit	0.094ng/ml	0.156-10ng/ml
	EH5293	Human Anti-Faricimab antibody ELISA Kit	0.469ng/ml	0.781-50ng/ml
Glargine	EH5242	Human Anti-Glargine antibody ELISA Kit	0.469ng/ml	0.781-50ng/ml
	EH5243	Human Anti-Glargine Neutralizing antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
Golimumab	EH5122	Human Golimumab ELISA Kit	0.375ng/ml	0.625-40ng/ml
	EH5109	Human Anti-Golimumab antibody ELISA Kit	0.188ng/ml	0.313-20ng/ml
Guselkumab	EU2694	Guselkumab ELISA Kit	0.469ng/ml	0.781-50ng/ml
	EU2695	Anti-Guselkumab Antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
Infliximab	EU3589	Infliximab ELISA Kit	0.188ng/ml	0.313-20ng/ml
	EU3590	Anti-Infliximab antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
Insulin Degludec	EH5324	Human Anti-Insulin Degludec antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
Integrin alpha V beta 6	EH5534	Human Anti-ITG $\alpha\beta 6$ Antibody ELISA Kit	0.469ng/ml	0.781-50ng/ml
Isatuximab	EH5327	Human Anti-Isatuximab antibody ELISA kit	0.469ng/ml	0.781-50ng/ml
Levilimab	EU2683	Levilimab (Anti-IL-6Ra) ELISA Kit	0.469ng/ml	0.781-50ng/ml

Drug	Cat.No	Product Name	Sensitivity	Range
Liraglutide	EU2691	LRT (Liraglutide) ELISA Kit	1.875ng/ml	3.125-200ng/ml
Natalizumab	EH5355	Human Anti-Natalizumab antibody ELISA Kit	0.94ng/ml	1.56-100ng/ml
Nivolumab	EH5349	Human Nivolumab ELISA Kit	0.469ng/ml	0.781-50ng/ml
	EH5348	Human Anti-Nivolumab antibody ELISA Kit	0.94ng/ml	1.56-100ng/ml
	EH5349-CM	Human Nivolumab ELISA Kit	0.938ng/ml	1.563-100ng/ml
Ocrelizumab	EU6141	Anti-Ocrelizumab antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
Omalizumab	EU2648	Omalizumab ELISA Kit	0.469ng/ml	0.781-50ng/ml
Panitumumab	EU2667	Panitumumab ELISA Kit	0.469ng/ml	0.781-50ng/ml
	EU2668	Anti-Panitumumab antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
Pembrolizumab	EU2647	Pembrolizumab ELISA Kit	0.938ng/ml	1.563-100ng/ml
	EU2671	Anti-Pembrolizumab antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
	EU2712	Anti-Pembrolizumab Neutralizing antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
Pertuzumab	EH5145	Human Pertuzumab ELISA Kit	0.469ng/ml	0.781-50ng/ml
	EH5146	Human Anti-Pertuzumab antibody ELISA kit	0.938ng/ml	1.563-100ng/ml
	EU2693	Anti-Pertuzumab Neutralizing Antibody ELISA kit	4.688ng/ml	7.813-500ng/ml
Ramucirumab	EH5310	Human Ramucirumab ELISA kit	0.938ng/ml	1.563-100ng/ml
Recombinant	EH5113	Recombinant Anti-Human CD28 mAb ELISA Kit	0.281ng/ml	0.469-30ng/ml
	EH5114	Recombinant Anti-Human CD3 ϵ mAb ELISA Kit	0.188ng/ml	0.313-20ng/ml
Rituximab	EH5295	Human Anti-Rituximab antibody ELISA Kit	0.469ng/ml	0.781-50ng/ml
	EU2616	Rituximab ELISA Kit	0.188ng/ml	0.313-20ng/ml
Romosozumab	EH5131	Human Romosozumab ELISA Kit	0.469ng/ml	0.781-50ng/ml
	EH5292	Human Anti-Romosozumab antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
Sacituzumab	EU6794	Sacituzumab govitecan (Anti-TROP2) ELISA Kit	0.469ng/ml	0.781-50ng/ml
Semaglutide	EH5240	Human Anti-Semaglutide antibody ELISA Kit	0.469ng/ml	0.781-50ng/ml
	EH5241	Human Anti-Semaglutide Neutralizing antibody ELISA Kit	0.469ng/ml	0.781-50ng/ml
	EH5258	Human Semaglutide Immunoglobulin E ELISA Kit	0.938ng/ml	1.563-100ng/ml
	EU2706	Semaglutide QuickTest ELISA Kit	56.25pg/ml	93.75-6000pg/ml
Teduglutide	EH5335	Human Anti-Teduglutide antibody ELISA kit	0.938ng/ml	1.563-100ng/ml
Tocilizumab	EH5315	Human Tocilizumab ELISA kit	0.469ng/ml	0.781-50ng/ml
Trastuzumab	EH5144	Human Trastuzumab ELISA Kit	0.469ng/ml	0.781-50ng/ml
	EH5252	Human Anti-Trastuzumab antibody ELISA Kit	0.469ng/ml	0.781-50ng/ml
	EH5357	Human Anti-T-DXd antibody ELISA Kit	0.94ng/ml	1.56-100ng/ml
Trastuzumab	EU2672	Trastuzumab Emtansine ELISA Kit	0.469ng/ml	0.781-50ng/ml
	EH5252-CM	Human Anti-Trastuzumab antibody ELISA Kit	6ng/ml	10-640ng/ml
Ustekinumab	EH5321	Human Anti-Ustekinumab antibody ELISA Kit	0.938ng/ml	1.563-100ng/ml
	EH5320	Human Ustekinumab ELISA kit	0.938ng/ml	1.563-100ng/ml
Vedolizumab	EH5322	Human Anti-Vedolizumab antibody ELISA Kit	0.469ng/ml	0.781-50ng/ml
VEGF	EH5142	Anti-Human VEGF (Humanized Antibody)ELISA Kit	0.188ng/ml	0.313-20ng/ml