

anti- NDRG3 antibody

Product Information


Catalog No.:	FNab09845
Size:	100 μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20 $^{\circ}$ C for 24 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	N-myc downstream-regulated gene 3 protein
Synonyms:	Protein NDRG3,N-myc downstream-regulated gene 3 protein
Calculated MW:	50 kDa
Uniprot ID :	Q9UGV2

Application

Specificity:	Human, Mouse; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500 - 1:2000
Image:	



mouse brain were subjected to SDS PAGE followed by western blot with FNab09845(anti- NDRG3 antibody) at dilution of 1:1000