

---

## FITC Anti-Human CD56 (B-A19)

### Product Information

Catalog No.:	FITC-30008
Conjugation:	FITC
Size:	50T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	B-A19
IsoType:	IgG1, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Synonyms:	Neural cell adhesion molecule 1, Leu-19, NKH1, NCAM-1, NCAM1, NCAM
Uniprot ID :	P51677

### Application

Reactivity:	Human
Tested Application:	FC
Recommended dilution: Volume per test:	5 $\mu$ L