
HRP-conjugated 6X His tag monoclonal antibody

Product Information

Catalog No.:	FNSA-0019
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was purified from antisera by Protein A+G purification
Conjugation:	HRP purified from horseradish root
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 3mg/ml BSA, 50% glycerol, pH7.6. - 20°C for 24 months , avoid repeated freezing and thawing.

Application

Tested Application:ELISA, WB

Suggested dilution: 1:1000-1:10000 for western blotting with ECL substrates

1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates

Note: This product is for research use only, not for diagnostic or therapeutic use.