

AP-conjugated Affinipure Donkey Anti-Human IgG (H+L)

Product Information

| | |
|---------------|---|
| Catalog No.: | FNSA-0026 |
| Size: | 100μl/500μl |
| Form: | liquid |
| Purification: | The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Conjugation: | Alkaline Phosphatase purified from calf intestine |
| Stabilizer: | 5mg/ml BSA (IgG-Free, Protease-Free) |
| Storage: | 0.01M Tris-HCl, 0.25M NaCl, 50% glycerol, pH8.0. -20°C for 24 months , avoid repeated freezing and thawing. |

Application

Tested Application:ELISA, WB

Suggested dilution: 1:1000-1:3000 for ELISA and Western blotting

Note: This product is for research use only, not for diagnostic or therapeutic use.