

---

## Alexa Fluor 488-Goat Anti-Rabbit IgG(H+L)

### Product Information

Catalog No.:	FNSA-0056
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was purified by immunoaffinity chromatography
Conjugation:	Alexa Fluor 488
Fluorophore:	Alexa Fluor 488 carboxylic acid, Amax=493nm; Emax=519nm
Host:	Goat
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.

### Application

Tested Application:	IF, FC
Recommended dilution:	1:100-1:500

**Note: This product is for research use only, not for diagnostic or therapeutic use.**

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone, Wuhan, Hubei, China (430206)

Tel: (0086) 027-87384275

Fax: (0086) 027-87800889

[www.fn-test.com](http://www.fn-test.com)