

---

## Alexa Fluor 594 – conjugated Donkey Anti-Mouse IgG(H+L)

### Product Information

Catalog No.:	FNSA-0062
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was purified by immunoaffinity chromatography
Conjugation:	Alexa Fluor 594
Fluorophore:	Alexa Fluor 594 carboxylic acid
Fluorophore/Protein Absorbance Ratio:	A591/A280=0.83
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.

### Application

Tested Application:	IF, FC
Recommended dilution:	1:100-1:500

**Note: This product is for research use only, not for diagnostic or therapeutic use.**