

Rhodamine (TRITC)–conjugated Goat Anti-Rabbit IgG (H+L)

Product Information

Catalog No.:	FNSA-0066
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was purified by immunoaffinity chromatography
Fluorophore:	Teramethyl Rhodamine Isothiocyanate(TRITC) A _{max} =550nm, E _{max} =570nm
Fluorophore/Protein Absorbance Ratio:	A ₅₅₀ /A ₂₈₀ =0.49
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.

Application

Tested Application:	IF, FC
Recommended dilution:	1:20-1:100

Note: This product is for research use only, not for diagnostic or therapeutic use.