

## **Rhodamine (TRITC)–conjugated Donkey Anti-Goat IgG (H+L)**

### **Product Information**

Catalog No.:	FNSA-0067
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was purified by immunoaffinity chromatography
Fluorophore:	Teramethyl Rhodamine Isothiocyanate(TRITC) A <sub>max</sub> =550nm, E <sub>max</sub> =570nm
Fluorophore/Protein Absorbance Ratio:	A <sub>550</sub> /A <sub>280</sub> =0.49
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.

### **Application**

Tested Application:	IF, FC
Recommended dilution:	1:20-1:100

**Note: This product is for research use only, not for diagnostic or therapeutic use.**