

## Cy3-conjugated Affinipure Goat Anti-Mouse IgG(H+L)

## **Product Information**

Catalog No.:	FNSA-0077
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was purified from antisera by immunoaffinity
	chromatography using antigens coupled to agarose beads.
Conjugation:	Purified from seaweed
Fluorophore:	Cy3 bis-NHS ester Amax=550nm,Emax=570nm
Fluorophore/Protein Absorbance Ratio: A552/A280=2.2	
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6.
	-20°C for 24 months, avoid repeated freezing and thawing.

## Application

Tested Application:	IF, FC
Recommended dilution:	1:20-1:100

## Note: This product is for research use only, not for diagnostic or therapeutic use.

Wuhan Fine Biotech Co., Ltd.B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone, Wuhan,<br/>Hubei, China (430206)Tel: (0086) 027-87384275Fax: (0086) 027-87800889www.fn-test.com