

Cy3-conjugated Affinipure Goat Anti-Rabbit IgG (H+L)

Product Information

Catalog No.: FNSA-0078

Size: $100 \mu l / 500 \mu l$

Form: liquid

Purification: The antibody was purified from antisera by immunoaffinity

chromatography using antigens coupled to agarose beads.

Conjugation: Purified from seaweed

Fluorophore: Cy3 bis-NHS ester Amax=550nm,Emax=570nm

Fluorophore/Protein Absorbance Ratio: A552/A280=2.2

Storage: 0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6.

-20°C for 24 months, avoid repeated freezing and thawing.

Application

Tested Application: IF, FC

Recommended dilution: 1:20-1:100

Note: This product is for research use only, not for diagnostic or therapeutic use.