
Rabbit anti Rhesus Monkey IgG

Product Information

Catalog No.:	FNSA-0092
Size:	100µg
Form:	liquid
Immunogen:	Rhesus Monkey IgG
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Specificity:	Rhesus Monkey
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 0.1% proclin300, pH7.4. 4°C for 24 months.

Note: This product is for research use only, not for diagnostic or therapeutic use.