

---

## Alexa Fluor 555-Goat Anti-Rabbit IgG(H+L)

### Product Information

Catalog No.:	FNSA-0095
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was purified by immunoaffinity chromatography
Conjugation:	Alexa Fluor 555
Fluorophore:	Alexa Fluor 555 carboxylic acid
Host:	Goat
Storage:	0.01M Sodium phosphate, 0.25M NaCl,0.02% sodium azide, 50% glycerol, pH7.6. -20°C, avoid repeated freezing and thawing. Store in dark.

### Application

Tested Application:	IF, FC
Recommended dilution:	1:100-1:1000

**Note: This product is for research use only, not for diagnostic or therapeutic use.**