
Normal Rabbit IgG Control

Product Information

| | |
|---------------|--|
| Catalog No.: | FNSA-0106 |
| Size: | 100µl/ 500µl |
| Form: | liquid |
| Source: | polyclonal rabbit IgG |
| Purification: | Protein A +G purified |
| Storage: | 0.01M Sodium phosphate, 0.25M NaCl, 0.1% proclin300, 50% glycerol, pH7.4. -20°C for 24 months , avoid repeated freezing and thawing. |

Application

Negative control for use in conjunction with other antibodies. Use at the same concentration as the detection antibody.

Note: This product is for research use only, not for diagnostic or therapeutic use.