
Rat IgG- Isotype Control

Product Information

Catalog No.:	FNSA-0134
Size:	1mg
Form:	Liquid
Clone:	polyclonal
Isotype:	IgG
Purity:	≥95% as determined by SDS-PAGE
Source:	Normal Rat IgG is obtained from rat serum that have not been Immunized.
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 0.1% proclin300, pH 7.4. 12 months from date of receipt at -20°C. Avoid freeze / thaw cycle.

Application

Negative control for use in conjunction with other antibodies.

Note: This product is for research use only, not for diagnostic or therapeutic use.