

Alexa Fluor 594 conjugated Affinipure Rabbit Anti-Rat IgG(H+L)

Product Information

Catalog No.:	FNSA-0135
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was purified by immunoaffinity chromatography
Conjugation:	Alexa Fluor 594
Fluorophore:	Alexa Fluor 594 carboxylic acid
Fluorophore/Protein Absorbance Ratio:	A591/A280=0.83
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.

Application

Tested Application:	IF, FC
Recommended dilution:	1:100-1:500

Note: This product is for research use only, not for diagnostic or therapeutic use.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone, Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com