

anti- 14-3-3 antibody

Product Information

Catalog No.:	FNab00001
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

14-3-3 theta(also known as YWHAQ, or 14-3-3 tau) is a member of 14-3-3 proteins which were the first phosphoserine/phosphothreonine-binding proteins to be discovered. 14-3-3 family members interact with a wide spectrum of proteins and possess diverse functions. Mammals express seven distinct 14-3-3 isoforms(gamma, epsilon, beta, zeta, sigma, theta, tau) that form multiple homo-and hetero-dimers. 14-3-3 proteins display the highest expression levels in the brain, and have been implicated in several neurodegenerative diseases, including Alzheimer's disease and amyotrophic lateral sclerosis. This antibody was raised against full-length 14-3-3 theta.

Immunogen information

Immunogen:	tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, theta polypeptide
Synonyms:	14 3 3, 14 3 3 protein T cell, 14 3 3 protein tau, 14 3 3 protein theta, 1C5, HS1, Protein HS1, YWHAQ
Observed MW:	28 kDa
Uniprot ID :	P27348

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com


Application

Reactivity: Human, Mouse, Rat


Tested Application: ELISA, IF, IP, IHC, WB

Recommended dilution: WB: 1:500-1:2000; IP: 1:200-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200

Image:



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using FNab00001(14-3-3 Antibody) at dilution of 1:50



Immunofluorescent analysis of HepG2 cells, using FNab00001(14-3-3 Antibody) at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



IP Result of anti 14-3-3(IP:FNab000001, 5μg;
Detection:FNab000001 1:500)with mouse lung
tissue lysate 4000μg.



HeLa cell were subjected to SDS PAGE followed
by western blot with FNab000001(14-3-3 Antibody)
at dilution of 1:1000