

anti- 5 Lipoxygenase antibody

Product Information

Catalog No.: FNab00007

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

ALOX5(arachidonate 5-lipoxygenase) is also named as LOG5, 5-LO and belongs to the lipoxygenase family. It catalyzes the first step in leukotriene biosynthesis, and plays a role in inflammatory processes. It is a 78 kDa iron-containing enzyme that translocates from the cytosol to associate with FLAP, and converts arachidonate in two steps into the unstable intermediate LTA4. It is essential for cysteinyl-leukotriene(cys-LT) production, critical mediators in asthma(PMID:12911785).

Immunogen information

Immunogen: arachidonate 5-lipoxygenase

Synonyms: 5 Lipoxygenase, 5 LO, 5 LOX, 5LPG, ALOX5, LOG5

Observed MW: 68-78 kDa
Uniprot ID: P09917

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, IF, WB

Recommended dilution: WB: 1:200-1:2000; IF: 1:20-1:100

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1



Image:



HL-60 cells were subjected to SDS PAGE followed by western blot with FNab00007(ALOX5 Antibody) at dilution of 1:600