

anti- A4GNT antibody

Product Information

Catalog No.:	FNab00015
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Necessary for the synthesis of type III mucin. Catalyzes the transfer of N-acetylglucosamine(GlcNAc) to core 2 branched O-glycans.

Immunogen information

Immunogen:	alpha-1,4-N-acetylglucosaminyltransferase
Synonyms:	Alpha4GnT,A4GCT
Observed MW:	39 kDa
UniprotID :	Q9UNA3

Application


Reactivity:	Human, Mouse, Rat	
Tested Application:	ELISA, WB, IHC	
Recommended dilution: WB: 1:500-1:5000; IHC: 1:20-1:200		
Image:		

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





-70kd

-55kd -40kd

-35kd

-25kd

-15kd

-10kd

Immunohistochemistry of paraffin-embedded human stomach slide using FNab00015(A4GNT Antibody) at dilution of 1:50

K-562 cells were subjected to SDS PAGE -100kd followed by western blot with FNab00015(A4GNT antibody) at dilution of 1:1000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Fax: (0086)027-87800889 www.fn-test.com Tel :(0086)027-87384275