

# anti- ABCA7 antibody

#### **Product Information**

Catalog No.: FNab00028

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

Plays a role in phagocytosis by macrophages of apoptotic cells. Binds APOA1 and may function in apolipoprotein-mediated phospholipid efflux from cells. May also mediate cholesterol efflux. May regulate cellular ceramide homeostasis during keratinocytes differentiation.

# Immunogen information

Immunogen: ATP-binding cassette, sub-family A(ABC1), member 7

Synonyms: ABCA7 Synonyms ABCA SSN, ABCA7, ABCX, Autoantigen SS N,

FLJ40025, Macrophage ABC transporter

Observed MW: 234 kDa UniprotID: Q8IZY2

### **Application**

Reactivity: Human, Mouse Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:5000

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com

1



# Image:



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00028(ABCA7 Antibody) at dilution of 1:1000

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com