

anti- ABCG4 antibody

Product Information

Catalog No.:	FNab00043
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May be involved in macrophage lipid homeostasis.

Immunogen information

Immunogen:	ATP-binding cassette, sub-family G(WHITE), member 4
Synonyms:	ABCG4, WHITE2
Observed MW:	60 kDa
Uniprot ID :	Q9H172


Application

Reactivity:	Human, Mouse, Rat	
Tested Application:	ELISA, WB, IHC	
Recommended dilution: WB: 1:200-1:2000; IHC: 1:20-1:200		
Image:		

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)





Immunohistochemistry of paraffin-embedded human liver cancer using FNab00043(ABCG4 antibody) at dilution of 1:50

U-937 cells were subjected to SDS PAGE followed by western blot with FNab00043(ABCG4 antibody) at dilution of 1:1000



Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com