

# anti- ABHD5 antibody

#### **Product Information**

Catalog No.: FNab00050

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

### **Background**

Lysophosphatidic acid acyltransferase which functions in phosphatidic acid biosynthesis(PubMed:18606822). May regulate the cellular storage of triacylglycerol through activation of the phospholipase PNPLA2(PubMed:16679289). Involved in keratinocyte differentiation(PubMed:18832586). Regulates lipid droplet fusion(By similarity).

### **Immunogen information**

Immunogen: abhydrolase domain containing 5

Synonyms: ABHD5, CDS, CGI58, IECN2, NCIE2

Observed MW: 39 kDa UniprotID: Q8WTS1

#### **Application**

Reactivity: Human, Mouse

Tested Application: ELISA, WB, IF, IHC

Recommended dilution: WB: 1:200-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>


1



## Image:



Immunohistochemistry of paraffin-embedded human skeletal muscle using FNab00050(ABHD5 antibody) at dilution of 1:50



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab00050(ABHD5 antibody) at dilution of 1:500