

anti- ABLIM2 antibody

Product Information

Catalog No.: FNab00055

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

May act as scaffold protein. May stimulate ABRA activity and ABRA-dependent SRF transcriptional activity.

Immunogen information

Immunogen: actin binding LIM protein family, member 2

Synonyms: AbLIM 2, ABLIM2, Actin binding LIM protein 2, DKFZp761F129,

FLJ26534, FLJ39684, KIAA1808

Observed MW: 68 kDa UniprotID: Q6H8Q1

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:500-1:5000; IHC: 1:20-1:200; IP: 1:500-1:5000

Image:

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)


Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com

1





Immunohistochemistry of paraffin-embedded human skeletal muscle slide using FNab00055(ABLIM2 Antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00055(ABLIM2 antibody) at dilution of 1:1000