

anti- ACAA2 antibody

Product Information

Catalog No.: FNab00059

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Abolishes BNIP3-mediated apoptosis and mitochondrial damage.

Immunogen information

Immunogen: acetyl-Coenzyme A acyltransferase 2

Synonyms: ACAA2, Beta ketothiolase, FLJ35992, DSAEC, FLJ95265

Observed MW: 42 kDa UniprotID: P42765

Application

Reactivity: Human, Mouse, Rat Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:200-1:2000; IHC: 1:20-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>





Immunohistochemistry of paraffin-embedded human lung cancer using FNab00059(ACAA2 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00059(ACAA2 antibody) at dilution of 1:500