

anti- ACAD10 antibody

Product Information

Catalog No.: FNab00060

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Acyl-CoA dehydrogenase only active with R-and S-2-methyl-C15-CoA.

Immunogen information

Immunogen: acyl-Coenzyme A dehydrogenase family, member 10

Synonyms: ACAD 10
Observed MW: 119 kDa
UniprotID: Q6JQN1

Application

Reactivity: Human, Mouse, Rat Tested Application: ELISA, WB, IHC, IF

Recommended dilution: WB: 1:200-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human pancreas using FNab00060(ACAD10 antibody) at dilution of 1:100



mouse spleen tissue were subjected to SDS PAGE followed by western blot with FNab00060(ACAD10 antibody) at dilution of 1:300