

# anti- ACADVL antibody

### **Product Information**

Catalog No.: FNab00068

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

## **Background**

Very long-chain acyl-CoA dehydrogenase (VLCAD) is one of four flavoproteins which catalyze the initial step of the mitochondrial b-oxidation spiral. It belongs to the acyl-CoA dehydrogenase family and is a homodimer of a 71-kDa polypeptide. The molecular mass of the nondenatured trypsinized VLCAD is 98 kDa, by gel filtration chromatography, indicating that it is a homodimer of the 48 kDa(tryptic digest) polypeptide(PMID:9461620). Defects in ACADVL are the cause of acyl-CoA dehydrogenase very long chain deficiency (ACADVLD). It has 2 isoforms(70 kDa and 68 kDa) produced by alternative splicing and a transit peptide.

### **Immunogen information**

Immunogen: acyl-Coenzyme A dehydrogenase, very long chain

Synonyms: ACAD6, ACADVL, LCACD, VLCAD

Observed MW: 48 kDa, 68-73 kDa

Uniprot ID: P49748

# **Application**

Reactivity: Human, Mouse, Rat

## 1

## Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)


Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>




Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:200-1:2000; IHC: 1:20-1:200

Image:



Immunohistochemistry of paraffin-embedded human liver cancer using FNab00068(ACADVL antibody) at dilution of 1:50



Various lysates were subjected to SDS PAGE followed by western blot with FNab00068(ACADVL antibody) at dilution of 1:1000

## Wuhan Fine Biotech Co., Ltd.