

## anti- ACAP2 antibody

### Product Information

Catalog No.:	FNab00070
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

GTPase-activating protein(GAP) for ADP ribosylation factor 6(ARF6).

### Immunogen information

Immunogen:	ArfGAP with coiled-coil, ankyrin repeat and PH domains 2
Synonyms:	ACAP2, Centaurin beta 2, CENTB2, CNT B2, KIAA0041
Observed MW:	88 kDa
UniprotID :	Q15057

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:2000; IP: 1:200-1:1000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab00070(ACAP2 antibody) at dilution of 1:300

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)