

anti- ACAP3 antibody

Product Information

Catalog No.: FNab00071

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

GTPase-activating protein for the ADP ribosylation factor family.

Immunogen information

Immunogen: ArfGAP with coiled-coil, ankyrin repeat and PH domains 3

Synonyms: ACAP3, Centaurin beta 5, CENTB5, Cnt b5, KIAA1716

Observed MW: 92 kDa UniprotID: Q96P50

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:5000; IF: 1:50-1:500

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel: (0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Jurkat cells were subjected to SDS PAGE followed by western blot with FNab00071(ACAP3 antibody) at dilution of 1:500