

anti- ACBD4 antibody

Product Information

Catalog No.: FNab00075

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Binds medium-and long-chain acyl-CoA esters and may function as an intracellular carrier of acyl-CoA esters.

Immunogen information

Immunogen: acyl-Coenzyme A binding domain containing 4

Synonyms: ACBD4, FLJ13322, FLJ90623, HMFT0700

Observed MW: 38 kDa
UniprotID: Q8NC06

Application

Reactivity: Human

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:5000; IF: 1:20-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





human kidney tissue were subjected to SDS PAGE followed by western blot with FNab00075(ACBD4 Antibody) at dilution of 1:1000

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>