

anti- ACK1 antibody

Product Information

Catalog No.: FNab00086

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

TNK2, also named as ACK1, belongs to the protein kinase superfamily and Tyr protein kinase family. It is a non-receptor tyrosine-protein and serine/threonine-protein kinase that is implicated in cell spreading and migration, cell survival, cell growth and proliferation. TNK2 transduces extracellular signals to cytosolic and nuclear effectors. It is a downstream effector of CDC42 which mediates CDC42-dependent cell migration via phosphorylation of BCAR1. TNK2 confers metastatic properties on cancer cells and promotes tumor growth by negatively regulating tumor suppressor such as WWOX and positively regulating pro-survival factors such as AKT1 and AR. The 60-70kd band is the isoform 2 of TNK2. Isoform 1 and 3 is about 120-125kd.

Immunogen information

Immunogen: tyrosine kinase, non-receptor, 2

Synonyms: TNK2, ACK1, ACK, FLJ45547, ACK 1

Observed MW: 70 kDa Uniprot ID: Q07912

Application

1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>




Reactivity: Human, Mouse, Rat


Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:200-1:2000; IHC: 1:20-1:200

Image:



Immunohistochemistry of paraffin-embedded human breast cancer using FNab00086(ACK1 antibody) at dilution of 1:50



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab00086(ACK1 antibody) at dilution of 1:500