

anti- ACOT6 antibody

Product Information

Catalog No.:	FNab00092
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Background	

Background

None

Immunogen information

Immunogen:	acyl-CoA thioesterase 6
Synonyms:	ACOT6, acyl CoA thioesterase 6, c14_5530, C14orf42
Observed MW:	23 kDa
Uniprot ID :	Q3I5F7

Application


Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB, IHC		
Recommended dilution: WB: 1:500-1:5000; IHC: 1:20-1:200			
Image:			

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com





Immunohistochemistry of paraffin-embedded human kidney slide using FNab00092(ACOT6 Antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab00092(ACOT6 antibody) at dilution of 1:1000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)