

anti- ACOX2 antibody

Product Information

Catalog No.:	FNab00096
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Oxidizes the CoA esters of the bile acid intermediates di-and tri-hydroxycholestanoic acids.

Immunogen information

Immunogen:	acyl-Coenzyme A oxidase 2, branched chain
Synonyms:	ACOX2, THCCox, THCA CoA oxidase, BRCOX, BRCACOX
Observed MW:	77 kDa
Uniprot ID :	Q99424

Application

Reactivity:	Human, Mouse, Rat	
Tested Application:	ELISA, IHC, IF	
Recommended dilution:IF: 1:10-1:100; IHC: 1:20-1:200		
Image:		


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com





Immunohistochemical of paraffin-embedded human liver using FNab00096(ACOX2 antibody) at dilution of 1:50

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com