

anti- ACOXL antibody

Product Information

Catalog No.:	FNab00097
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	acyl-Coenzyme A oxidase-like
Synonyms:	ACOXL, Acyl CoA oxidase like protein, acyl Coenzyme A oxidase like, FLJ11042
Observed MW:	32 kDa
Uniprot ID :	Q9NUZ1

Application

Reactivity:	Human
Tested Application:	ELISA, IP, IHC
Recommended dilution:	IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human spleen tissue slide using FNab00097(ACOXL Antibody) at dilution of 1:50



IP Result of anti-ACOXL (IP: FNab00097, 3ug; Detection: FNab00097 1:500) with K-562 cells lysate 2000ug.