

# anti- ACTL7B antibody

#### **Product Information**

Catalog No.:	FNab00120
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Background	

None

## **Immunogen information**

Immunogen:	actin-like 7B
Synonyms:	Actin like 7 beta, actin like 7B, Actin like protein 7B, ACTL7B
Observed MW:	45 kDa
UniprotID :	Q9Y614

## Application


Reactivity:	Human, Mouse, Rat			
Tested Application:	ELISA, WB, IHC			
Recommended dilution: WB: 1:500-1:5000; IHC: 1:20-1:200				
Image:				

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





-35kd

-25kd

Immunohistochemistry of paraffin-embedded human prostate cancer using FNab00120(ACTL7B antibody) at dilution of 1:100

human testis tissue were subjected to SDS PAGE -100kd -70kd 1:300 -55kd 40kd

followed by western blot with FNab00120(ACTL7B antibody) at dilution of

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com