

anti- ACTR10 antibody

Product Information

| Catalog No.: | FNab00125 |
|---------------|---|
| Cutulog 110 | 11/4000125 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | \geq 95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |
| Background | |

Background

None

Immunogen information

| Immunogen: | actin-related protein 10 homolog(S. cerevisiae) |
|--------------|---|
| Synonyms: | Actin related protein 10, Actin related protein 11, ACTR10, ACTR11, ARP11, HARP11 |
| Observed MW: | 46 kDa |
| UniprotID : | Q9NZ32 |

Application


| Reactivity: | Human, Mouse, Rat |
|-----------------------|---|
| Tested Application: | ELISA, WB, IHC, IF |
| Recommended dilution: | WB: 1:200-1:2000; IHC: 1:20-1:200; IP: 1:200-1:1000; IF: 1:10-1:100 |
| Image: | |

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human hepatocirrhosis using FNab00125(ACTR10 antibody) at dilution of 1:100

mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00125(ACTR10 Antibody) at dilution of 1:600



Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>