

# anti- ACTR3B antibody

#### **Product Information**

Catalog No.: FNab00127

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

### **Background**

Plays a role in the organization of the actin cytoskeleton. May function as ATP-binding component of the Arp2/3 complex which is involved in regulation of actin polymerization and together with an activating nucleation-promoting factor(NPF) mediates the formation of branched actin networks. May decrease the metastatic potential of tumors.

### **Immunogen information**

Immunogen: ARP3 actin-related protein 3 homolog B(yeast)

Synonyms: ACTR3B, Actin like protein 3B, ARP3 beta, ARP4, ARP11

Observed MW: 48 kDa UniprotID: Q9P1U1

#### **Application**

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:5000

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com

1



## Image:



human brain tissue were subjected to SDS PAGE followed by western blot with FNab00127(ACTR3B antibody) at dilution of 1:800

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>