

anti- ADAMDEC1 antibody

Product Information

Catalog No.: FNab00145

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

May play an important role in the control of the immune response and during pregnancy.

Immunogen information

Immunogen: ADAM-like, decysin 1

Synonyms: ADAM DEC1, ADAM like protein decysin 1, ADAM like, decysin 1,

ADAMDEC1, FLJ79219, M12.219

Observed MW: 53 kDa, 44 kDa

UniprotID: O15204

Application

Reactivity: Human

Tested Application: ELISA, WB, IHC

Recommended dilution: WB: 1:200-1:1000; IHC: 1:20-1:200; IP: 1:50-1:500

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab00145(ADAMDEC1 Antibody) at dilution of 1:50



MCF-7 cells were subjected to SDS PAGE followed by western blot with FNab00145(ADAMDEC1 Antibody) at dilution of 1:300