

anti- ADAMTS12 antibody

Product Information

Catalog No.: FNab00146

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Metalloprotease that may play a role in the degradation of COMP. Cleaves also alpha-2 macroglobulin and aggregan. Has anti-tumorigenic properties.

Immunogen information

Immunogen: ADAM metallopeptidase with thrombospondin type 1 motif, 12

Synonyms: ADAM TS 12, ADAM TS12, ADAMTS 12, ADAMTS12, PRO4389

Observed MW: 170 kDa UniprotID: P58397

Application

Reactivity: Human, Mouse

Tested Application: ELISA, WB, IHC, IF

Recommended dilution: WB: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200

Image:


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human stomach cancer tissue slide using FNab00146(ADAMTS12 Antibody) at dilution of 1:50



mouse placenta tissue were subjected to SDS PAGE followed by western blot with FNab00146(ADAMTS12 Antibody) at dilution of 1:300

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)