

# anti- ADAT2 antibody

### **Product Information**

Catalog No.:	FNab00153
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

#### Background

Probably participates in deamination of adenosine-34 to inosine in many tRNAs.

### **Immunogen information**

Immunogen:	adenosine deaminase, tRNA-specific 2, TAD2 homolog(S. cerevisiae)
Synonyms:	ADAT2, DEADC1, dJ20N2, dJ20N2.1, DKFZp686L1118, FLJ44213, TAD2
Observed MW:	22 kDa
UniprotID :	Q7Z6V5

## Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution	:WB: 1:200-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human hepatocirrhosis using FNab00153(ADAT2 antibody) at dilution of 1:50



A431 cells were subjected to SDS PAGE followed by western blot with FNab00153(ADAT2 antibody) at dilution of 1:400

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>