

# anti- ADCK5 antibody

#### **Product Information**

Catalog No.: FNab00154

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

### **Background**

The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate(Ser, Thr or Tyr).

#### **Immunogen information**

Immunogen: aarF domain containing kinase 5

Synonyms: ADCK5, FLJ35454

Observed MW: 66 kDa
UniprotID: Q3MIX3

## **Application**

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Recommended dilution: WB: 1:200-1:2000

Image:


#### 1

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab00154(ADCK5 antibody) at dilution of 1:800

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com