

anti- ADCY8 antibody

Product Information

Catalog No.: FNab00158

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

This is a membrane-bound, calcium-stimulable adenylyl cyclase. May be involved in learning, in memory and in drug dependence(By similarity).

Immunogen information

Immunogen: adenylate cyclase 8(brain)

Synonyms: AC8, ADCY3, ADCY8, adenylate cyclase 8(brain), Adenylate cyclase

type 8, Adenylate cyclase type VIII, Adenylyl cyclase 8, ATP

pyrophosphate lyase 8, HBAC1

Observed MW: 140 kDa UniprotID: P40145

Application

Reactivity: Human, Mouse Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:5000

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com

1



Image:



Neuro-2a cells were subjected to SDS PAGE followed by western blot with FNab00158(ADCY8 Antibody) at dilution of 1:1000

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>