

anti- ADH1B antibody

Product Information

Catalog No.:	FNab00162
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	\geq 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Background	

None

Immunogen information

Immunogen:	alcohol dehydrogenase 1B(class I), beta polypeptide
Synonyms:	ADH1B, ADH2, Alcohol dehydrogenase 1B, DKFZp686C06125
Observed MW:	40 kDa
UniprotID :	P00325

Application


Reactivity:	Human, Mouse, Rat			
Tested Application:	ELISA, WB			
Recommended dilution:WB: 1:500-1:5000				
Image:				

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com