

anti- ADH7 antibody

Product Information

Catalog No.: FNab00166

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Could function in retinol oxidation for the synthesis of retinoic acid, a hormone important for cellular differentiation. Medium-chain(octanol) and aromatic(m-nitrobenzaldehyde) compounds are the best substrates. Ethanol is not a good substrate but at the high ethanol concentrations reached in the digestive tract, it plays a role in the ethanol oxidation and contributes to the first pass ethanol metabolism.

Immunogen information

Immunogen: alcohol dehydrogenase 7(class IV), mu or sigma polypeptide

Synonyms: ADH 4, ADH7, Gastric alcohol dehydrogenase, Retinol dehydrogenase

Observed MW: 41 kDa, observed 45-47 kDa

UniprotID: P40394

Application

Reactivity: Human

Tested Application: ELISA, WB, IF, IHC

Recommended dilution: WB: 1:500-1:5000; IF: 1:20-1:200; IHC: 1:20-1:200

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com


1



Image:



Immunohistochemistry of paraffin-embedded human liver slide using FNab00166(ADH7 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab00166(ADH7 antibody) at dilution of 1:1000