

anti- AF-6 antibody

Product Information

Catalog No.:	FNab00190	
Size:	100µg	
Form:	liquid	
Purification:	Immunogen affinity purified	
Purity:	\geq 95% as determined by SDS-PAGE	
Host:	Rabbit	
Clonality:	polyclonal	
Clone ID:	None	
IsoType:	IgG	
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)	

Background

MLLT4, also named as AF6 and Afadin, belongs to an adhesion system, probably together with the E-cadherin-catenin system, which plays a role in the organization of homotypic, interneuronal and heterotypic cell-cell adherens junctions(AJs). It is a nNectin-and actin-filament-binding protein that connects nectin to the actin cytoskeleton. The antibody is specific to MLLT4.

Immunogen information

Immunogen:	myeloid/lymphoid or mixed-lineage leukemia(trithorax homolog, Drosophila); translocated to, 4
Synonyms:	AF 6, AF6, AFADIN, Afadin variant, MLLT4
Observed MW:	207 kDa
Uniprot ID :	P55196

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889

www.fn-test.com



Recommended dilution:WB: 1:500-1:2000 Image:



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab00190(MLLT4 antibody) at dilution of 1:500

2

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

www.fn-test.com