

anti- AF9 antibody

Product Information

Catalog No.: FNab00191

Size: 100μg Form: liquid

Purification: Immunogen affinity purified

Purity: ≥95% as determined by SDS-PAGE

Host: Rabbit

Clonality: polyclonal

Clone ID: None IsoType: IgG

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen: myeloid/lymphoid or mixed-lineage leukemia(trithorax homolog,

Drosophila); translocated to, 3

Synonyms: AF9, FLJ2035, MLLT3, Protein AF 9, YEATS3

Observed MW: 64 kDa UniprotID: P42568

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB

Recommended dilution: WB: 1:50-1:1000

Image:


1

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Jurkat cells were subjected to SDS PAGE followed by western blot with FNab00191(MLLT3 antibody) at dilution of 1:600