

## anti- AFAP1L1 antibody

### Product Information

Catalog No.:	FNab00193
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

May be involved in podosome and invadosome formation.

### Immunogen information

Immunogen:	actin filament associated protein 1-like 1
Synonyms:	AFAP1 like protein 1, AFAP1L1, FLJ36748
Observed MW:	86 kDa
Uniprot ID :	Q8TED9

### Application

Reactivity:	Human, Monkey
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:2000
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



COS-7 cells were subjected to SDS PAGE followed by western blot with FNab00193(AFAP1L1 Antibody) at dilution of 1:600